Questa documentazione deve essere fornita all'utilizzatore del complesso tubo-guaina  
The contents of this documentation must be transmitted to the user of the tube-assembly  
Le contenu de cette documentation doit être transmis à l'utilisateur de la gaine équipée

Documentazione N° Documentation N° N° de Documentation	Revisione Edition Version	Data di edizione Date of release Date de l'édition	Testo originale Original text Texte original
INFO_TG	G	19.02.15	italiano / italian / italien

## 0. Sommario - Table of contents - Table des matières

0.	Sommario - Table of contents - Table des matières .....	2
1.	Destinazione d'uso - Use destination - Destination d'emploi.....	3
2.	Descrizione del funzionamento - Description of operation - Description du fonctionnement .....	3
3.	Aspetti generali - Common aspects - Aspects généraux.....	4
4.	Interfacciamento e compatibilità - Interfacing and compatibility - Interfacement et compatibilité .....	6
4.1	Collimatore - Collimator - Collimateur .....	6
4.2	Cavi - Cables - Câbles .....	6
4.3	Generatore e tavolo di comando - Generator and control board - Générateur et pupitre.....	7
4.4	Starter - Starter - Démarrleur.....	8
4.5	Supporti meccanici - Mechanical supports - Supports mécaniques .....	9
5.	Trasporto - Transportation - Transport.....	9
6.	Immagazzinamento - Storage - Stockage.....	9
7.	Installazione - Installation - Installation .....	10
7.1	Montaggio meccanico - Mechanical assembling - Montage mécanique .....	10
7.2	Cablaggio, cavi AT - Wiring, high voltage cables - Câblage, câbles haute tension.....	10
7.3	Verifiche: valore di preaccensione, accensione filamenti, rotazione - Preliminary checks: filament preheating, filament heating, rotation - Contrôles: préchauffage filaments, chauffage filaments, stator .....	11
7.4	Stabilità AT, formazione sull'impianto - High voltage stability, tube seasoning on the installation - Stabilité en haute tension, formation sur l'installation .....	12
7.5	Taratura - Calibration - Calibration.....	12
8.	Uso - Use - Emploi .....	12
8.1	Preriscaldamento giornaliero - Daily warm-up - Préchauffage journalier.....	13
8.2	Pause di raffreddamento - Cooling delays - Pauses de refroidissement.....	13
8.3	Limitazione del numero di lanci - Limitation of rotor acceleration cycles - Limitation des démarrages.....	14
8.4	Precauzioni in caso di versamento di liquidi - Precautions in case of accidental fluid spills - Précautions en cas de renversement de liquides.....	14
8.5	Pulizia esterna - External cleaning - Nettoyage extérieur .....	14
8.6	Formazione dopo periodi di non funzionamento - Seasoning after inoperative periods - Formation après des périodes de non fonctionnement .....	14
8.7	Superfici in temperatura - Hot surfaces - Surfaces en température .....	14
8.8	Bloccaggio della cravatta - Locking of trunnion rings - Blocage du collier de fixation.....	15
8.9	Rilevamento di guasti - Failure detection - Détection des pannes .....	15
9.	Manutenzione - Maintenance - Entretien.....	15
9.1	Verifica delle teste cavo - High voltage cables terminals - Embases haute tension .....	15
9.2	Verifica della stabilità in tensione - High voltage stability - Stabilité en tension .....	16
9.3	Verifica e regolazione della taratura - mA calibration - Verification et réglage de la calibration .....	16
9.4	Pulizia del ventilatore - Fan cleaning - Nettoyage du ventilateur .....	16
9.5	Verifica dei collegamenti BT e del conduttore di terra - Low voltage wiring and grounding - Câblages basse tension et conducteur de terre .....	16
9.6	Verifica delle sicurezze termiche - Thermal safety devices - Sécurités thermiques .....	16
9.7	Verifica delle masse sospese - Suspended masses - Masses suspendues .....	16
10.	Rimozione e smaltimento - Removal and disposal - Enlèvement et élimination .....	16
10.1	Complessi tubo guaina per diagnostica - Diagnostic tube assemblies - Gaine équipée pour diagnostic .....	16
10.2	Complessi tubo guaina per mammografia - Mammography tube assemblies - Gaine équipée pour mammographie .....	17
11.	Rischi residui - Residual risks - Risques résiduels .....	17
11.1	Radiazioni - Radiations - Radiations .....	17
11.1.1	Radiazione utile - Useful radiation - Radiation utile .....	17
11.1.2	Filtrazione insufficiente - Insufficient filtration - Filtration insuffisante .....	18
11.1.3	Radiazione di fuga - Leakage radiation - Radiation de fuite .....	18
11.2	Caduta della guaina - Tube assembly drop - Chute de la gaine .....	18
11.3	Alte tensioni - High voltages - Hautes tensions .....	18
11.4	Scariche verso massa delle parti in alta tensione - High voltage arcs to ground - Claquages vers la masse des parties en haute tension .....	19
11.5	Temperatura delle superfici esterne - External surfaces temperature - Température des surfaces externes .....	19
11.6	Esplosione della guaina - Housing explosion- Explosion de la gaine .....	20
11.7	Dimensioni errate dei fuochi - Incorrect focal spot size - Dimensions incorrectes des foyers .....	20
11.8	Rischio accessibilità guaina XK1016T – Accessibility hazard of XK1016T tube assembly – Risque accessibilité de la gaine équipée XK1016T .....	20
11.9	Rischio alimentazioni errate XK1016T – Incorrect power supply hazard of XK1016T – Risque alimentation incorrecte de la gaine équipée XK1016T.....	20
12.	Formazione sull'impianto - On site tube seasoning - Formation du tube sur l'installation .....	21

12.1	Riscaldamento - Warm-up - Préchauffage .....	21
12.2	Salita in kV - kV step up - Montée en kV.....	22
12.3	Scopia - Fluoroscopy - Scopie .....	22
12.4	Grafia - Radiography - Graphie .....	22
12.5	Note - Notes - Notes.....	22
12.6	Tabelle: Fasi di formazione - Seasoning phases – Phase de formation .....	24
13.	Scambiatore di calore locale - Local heat exchanger - Echangeur de chaleur local .....	30
13.1	Descrizione - Description - Description.....	30
13.2	Installazione - Installation - Installation.....	30
13.3	Manutenzione - Maintenance - Maintenance .....	31

## 1. Destinazione d'uso - Use destination - Destination d'emploi

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X-ray tubes assembly manufactured by IAE are intended for use in medical diagnostic, in association with the equipment designed for the same purpose.

They are not designed for use in sterile environment, exposed to atmospheric agents or in presence of deflating atmospheres.

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The X-ray tube is a vacuum tube operating from an electrical viewpoint like a diode in conditions of filament saturation.

When high voltage is applied between anode and cathode, the X-ray tube emits X radiation in all the half space above the tangent plane to the anode surface.

For this reason the tube must be used within a shielded envelope, the housing, which limits all the undesired part of X radiation, allowing exit to the useful X-ray beam only, through the housing output port.

The useful beam must be further limited by a fixed or adjustable collimator, which is not a part of the housing and must ensure, together with the housing, the continuity of leakage radiation shielding.

The housing also ensures protection from the high voltage which is applied to anode and cathode by the high voltage cables coming from an elevator transformer with rectifying bridge, usually called generator.

X-ray production requires the use of high energy, most of which is converted to heat dissipated into the tube elements.

The energy sent to the tube initially concentrates on the focal spot (milliseconds), then it distributes over the focal spot track (fractions of a sec-

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ond), then accumulates in the anode mass (seconds to minutes), later is transmitted by heat radiation to the insulating oil which is contained in the tube housing and to the housing metal walls, and from these to the ambient air by natural convection or forced ventilation (fractions of hours).

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These heat exchanges bring the tube elements to extreme temperatures, which must be kept within given limits to avoid destruction of the X-ray tube.

Focal spot and focal spot track temperatures are controlled by power ( $kV \times mA$ ) and exposure time (sec), which are supplied by the generator under the limitation of maximum load factors, as calculated from the rating charts.

Anode temperature is governed by a correct exposure to cooling delays ratio, as a function of the exposure energy.

The control of cooling delays must be ensured by the equipment resident software, by accumulated heat simulation circuitry or, in absence of these, by a careful and conservative programming of working schedules made by the user, on the basis of anode heating and cooling curves.

Housing temperature is also controlled by work to cooling delays ratio, by the same means described for the anode. In this case the duty cycles are calculated over longer intervals, of the order of half a working day, on the basis of housing heating and cooling curves. Energy computation must include the energy which is sent to the stator for rotor acceleration, which can be of the order of 10 kJ per cycle in high speed operation.

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Up to three safety devices are present in the tube unit for the housing temperature only: an external thermal switch, an internal one not accessible to the user and a microswitch on oil expansion bellows.

These devices are not to be considered as a substitute to automatic or manual programming of cooling delays, but as last safety levels never to be reached in normal operation.

**It is important to point out that no possible safety device exist in the tube unit for the other limiting temperatures, concerning focal spot, focal spot track and anode. For these, patient safety and tube integrity must be ensured by a correct control of energy cycles from the generator.**

Anode current as a function of filament current is given by the emission curves.

For construction reasons the tolerance on the filament current value which gives a determined value of anode current is comparatively wide (approx. 3% of filament current). For this reason the emission curves can be used as an initial operation setup, but cannot replace an individual tube calibration at each installation.

Pre-programmed systems containing information on filament emission must also allow for the same range of filament calibration.

The useful beam is filtered by the tube unit structural materials (glass, oil, plastic materials). However this inherent filtration is not sufficient to limit the low energy fraction of X radiation, the so called soft X-rays, which are harmful to the patient and do not contribute to image creation. This happens even more with Beryllium window tubes, typically mam-

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mography tubes, which emit a high quantity of soft X-rays.

For this reason an additional filtration must be placed on the useful beam, ensuring the compliance to the pertinent regulations while keeping account of all the interposed materials (collimator, compressor etc.)

A part of added filtration is normally present on the tube unit, but the obtainment of a sufficient total filtration is a task of the installer.

X-rays are anyway harmful to the patient. Since the diagnostic benefit to harm ratio can only be evaluated by a physician, X-ray equipment are to be operated exclusively by competent personnel under the guide of a radiologist.

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The X-ray tube assembly operates in conjunction with the following components of the radiological equipment, whose compatibility must be checked.

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Beam limiting devices must interface mechanically with the tube unit so that the continuity of leakage radiation shielding is ensured, and must have a leakage radiation low enough to comply IEC 601-1-3 or other local regulations, when added to the tube unit leakage radiation.

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For general diagnostic tube assemblies, suitable cables are shielded high voltage cables with three pins Federal cable heads, complying to the corresponding regulations and having a maximum length of 24 m.

For mammography tube assemblies, the suitable cable heads are the Alden type.

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Both high frequency and mains frequency generators can be used.

The X-ray tube assembly is a passive component not being able to limit the energy sent by the generator. For this reason, the generator - control board system must comply to the pertinent regulations, requiring that in any failure condition an uncontrolled energy supply to the tube is prevented as this eventuality could cause the tube assembly explosion with severe danger to the patient and the user.

The system must provide for safety circuits preventing high voltage applications to the tube insert in the event of activation of thermal safety devices.

**These circuits must be galvanically separated from mains, by an insulation complying with IEC 601.1 standard.**

Up to three safety devices, depending on the tube housing type, protect the tube unit from housing and oil overheating:

- a thermal switch on housing external surface, accessible to the installer.
- a limit switch on oil expansion bellows, accessible to the installer.
- an internal thermal switch, in series to the stator common wire. In case of non connection or failure of the first two, it causes the intervention of the stator supply safety devices or the tube failure.

These devices are to be intended anyway as a last safety level, never to be reached during the tube assembly normal use.

The limitation of the tube unit working temperature must thus be obtained by the following means:

- a suitable planning of exposure techniques and the cooling delays, in case of manual installations.
- a correct programming of automatic circuits and software, in case of automatic equipment.

Also, for physical reasons the tube as-

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sembly safety devices cannot be effective:

- in case of anode overheating (this occurrence must be prevented by cooling delays planning or by automatic systems, as discussed above)
- in case that a failure of the equipment makes it not responding to the intervention of tube housing safety devices (this event must be made impossible by an equipment design complying to the equipment regulations).

When designing a new equipment or replacing an existing tube assembly with a new or different make, it must be checked that the filament heating circuit can supply the filament voltage and current as required by filament characteristics and emission curves, including the calibration range required by filament current tolerances indicated on the same curves.

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The stator circuit must be able to supply the voltages and currents as detailed in the stator data sheet, and detect short or open circuits in the stator circuits, inhibiting in this case the high voltage to the tube.

**Stator power supply circuit must be electrically isolated from mains, by an insulation complying with IEC 601.1 standard.**

The stator boost time must be limited to the time required to attain the nominal speed, then be decreased to the stand-by value, with the purpose of limiting useless heating of the stator and tube unit.

For the same reason the number of boost cycles must be limited to the value mentioned in the data sheet, preferably by means of timing circuits. The danger of stator heating because of repeated boost cycles is particularly present when calibrating the equipment. The stator may be permanently damaged by a careless repeated boosting during calibration, creating the conditions for potential malfunctioning during normal operation. It is thus useful that at the tube design and installation instructions level the possibility of switching to continuous rotation during calibration be

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considered. This is also true for use in normal operation at exposure rates above one exposure per minute.

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Tube assembly support elements must comply to the safety requirements of IEC 601-1 standard, referring to suspended weights. The screws must be of class 12.9, and long enough to couple to all the available threads of tapped holes.

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#### 5. Trasporto - Transportation - Transport

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Tube assemblies, both new and returned, must always be shipped in the original packing, after checking that the packing is in good conditions.

While the housing structure is made by sturdy metal parts, the fragile part is the glass insert inside.

It is thus important that during all the transportation phases the tube unit be not subject to shocks and be kept in vertical position, the situation in which the glass part is less sensitive to shocks.

Should the packing be open for inspection and then undergo further transportation, care must be taken that all the packing parts are restored in their original position.

Upon receipt of the tube, the shock and tilt indicators should be checked for activation, and the subsequent actions taken according to the indicators instructions.

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Tube assemblies are to be stored in a non oxidizing, non corrosive, dry environment, within the temperature limits indicated in the data sheets.

Should the tubes be kept in stock for more than a year, the factory or the local representative should be contacted in order to repeat the tube test.

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According to the tube types, the tube unit can be fastened to the equipment by four M6 tapped holes in the window plane, four M6 or M8 tapped holes in the built-in mounting plate or by trunnion rings supplied as an add on.

Assembling must be done under the consideration that the tube assembly is a weight which may be suspended above the patient.

As discussed above, the mounting screws must be of class 12.9, and long enough to engage in all available threads. Should one of the threaded couplings appear damaged, or show signs of yielding when being tightened, the installation must be suspended until the problem is removed or the damaged parts are replaced.

Check that the parts connecting the shielding of the tube unit and the collimator must be integer and correctly in place.

Check that the total permanent filtration on the useful beam (tube unit and collimator inherent filtration + added filtration) comply with IEC standard 601.1.3 or the pertinent regulations applied in the country of installation.

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Prior than inserting the high voltage cables into the cable sockets, both these and the cable terminals must be thoroughly cleaned. Carbon traces from previous arcings, if any, must be completely removed. Cables must be replaced if they are darkened or swollen near the cable terminals, or if these show any crack. Check the integrity of the terminal pins, and restore their spring action if needed.

For each cable, build up a thick layer of silicon grease on the terminal extremity, and coat the lateral surface completely. Insert the cable termination into the socket, find the locating key, insert the cable nut into the fist

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threads and tighten slowly. The purpose is to draw the grease from the terminal extremity up to the outside, pushing out all the air remaining between the cable termination and the socket. Complete tightening at intervals during calibration, to compensate for grease setting with heating.

Always connect both thermal safety devices existing on the outside of tube unit. Both are a necessary safety element for the patient and the user.

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Check by the measurement of the filament transformer primary current and coil ratio, that the filament preheating current is between 2 A and 3 A. Too high a preheating current makes the filament evaporate, metallizes the glass and shortens the tube life.

Check that in radiography preparation the large and small filaments are selected correctly, to prevent focal track overloading in case of wrong focal spot selection.

WARNING to perform this check the high voltage to the tube must be inhibited, to avoid the risk of being irradiated from a short distance.

Make a shortened rotation boost to check the correct rotation direction, i.e. anode rotating CCW when looking at the anode from the cathode side.

Check that the rotor full speed is attained by the normal boost cycle. In RTM tubes the rotor speed can be measured by the frequency of the rotor feedback voltage in the stator windings. To perform the measurement, the starter has to be disconnected from the stator, allowing the rotor to rotate by inertia. In these conditions a voltage of approx. 400 mV, at the rotation frequency, is present between the stator common wire and each the other two wires. (RPM =  $60 \times f$  (Hz)).

In the other tubes the speed can be determined by a vibrometer, a stroboscope or a reed frequency meter. A non quantitative indication can be obtained by the rotation noise. If the boost cycle does not attain full rotor

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speed, the noise frequency will increase during rotation stand-by.

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If the last tube test was earlier than 2 months, perform the seasoning procedure on the installation according to the supplied instructions.

The same procedure must be performed if the tube has a tendency to make arcs.

In any case perform the warm up section of this procedure.

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Perform the calibration according to the equipment instructions.

To prevent tube unit overheating and permanent stator damage, the number of boost cycles must be limited to the value indicated in the data sheets. To avoid an excessive lengthening of the calibration procedure, it may be useful to switch the stator power supply to the stand-by voltage.

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In normal operation the user has never access to the tube assembly, because the operation is always controlled by the equipment control board. For this reason, specific instructions for the tube assembly cannot be given. Reference should be made to the equipment instructions.

Below some general purpose indications are given, being relevant or safety of use or tube life.

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Applying high loads to the tube with cold anode can prematurely wear the focal track, because the anode material has a fragile behaviour at room temperature, and thus its surface is more easily cracked under the thermal cycles induced by the load application.

While feasibility or emergency reasons are often predominant in diagnostic use, it is however important for tube life to perform a warm-up, as it is usual in computed tomography installations, if the tube has been left inoperative for more than 3 hours and the tube is operated at the highest powers. If the warm-up procedure is not included in the equipment software, the necessary data can be obtained by the first part of tube seasoning procedure here included.

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While the control board normally limits the peak radiographic power (maximum loads) this not always happens for the anode and housing thermal limits (accumulated energy and continuous power dissipation). However these are extremely important for a safe use of the equipment.

In absence of an equipment resident software having this purpose, the daily working schedule must be determined, with the help of the installation engineer if needed, so that:

- the energy of one examination is lower than the maximum anode heat content
- the average power of consecutive exams (work + delays) is lower than the maximum anode dissipation
- the average power of half a working day is lower than the maximum housing dissipation.

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Normally every time the first step of the exposure switch is depressed (exposure preparation), a rotor acceleration cycle is performed and the filament is overheated.

To prevent stator and housing overheating, the number of acceleration cycles must be limited (2 per minute at low speed, 1 per minute at high speed)

To extend filament life, avoid staying in the exposure preparation phase longer than needed.

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Live electrical circuits are located under the housing non-waterproof end covers.

In case of accidental fluid spills, cut power to the equipment immediately, then remove the covers, wipe off any moisture and replace the covers prior than switching the power on again.

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The housing external surface may be cleaned by a cloth damped with a water based detergent solution.  
**Power must be cut off during cleaning.**

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If the equipment has remained inoperative:

- for more than a month, perform the warm-up cycle
- for more than 6 months, perform the seasoning procedure.

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Because of the need of dissipating high energies towards the surrounding environment, the housing external surfaces can be heated up to 80 °C. The patient and the user must be prevented from coming in contact with the housing.

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When mounted on trunnion rings, the housing can be rotated for orientation of the X-ray beam and then be locked in position by a brake.

The brake must be tightened carefully to prevent an accidental rotation of the housing - collimator system which may hit the patient.

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Normally tube assembly failures (high voltage arcs, filament cut off, stator failures) are detected by the equipment safety devices.

However it is important to act at a prevention level every time any abnormal noise coming from the tube assembly or other parts of the high voltage circuit (cables, generator) may give suspect of a state of failure. Under these circumstances **suspend exposure immediately, cut off power to the equipment and remove the patient**, then contact the service engineer and follow his instruction.

Also perform an immediate action, either directly or through the service engineer, if any suspended mass (housing or collimator) shows an unstable fastening.

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In the occasion of ordinary equipment maintenance, the following checks are required by the tube unit.

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Check for presence of enough grease or oil in the cable receptacles, and absence of arc paths.

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By means of an oscilloscope and a high voltage probe, either invasive or not, that no part of the high voltage circuit, including the tube assembly, show a tendency to arcs.

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Clean the parts crossed by cooling air flow if present.

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Check that the wiring is efficient and that the grounding connection complies with the pertinent regulations.

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Check that all safety devices are connected and that opening the safety circuit prevents high voltage from being applied to the tube assembly.

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Check the screws tightening and the mechanical supports efficiency.

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The tube assembly contains the following materials, whose disposal must be done in respect of the local regulations:

- PCB free insulating oil
- Tungsten, molybdenum, rhenium, nickel, copper, silver
- Glass.
- Iron, aluminium and copper alloys.

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- Lead, lead oxide, epoxy resin, insulating materials, synthetic rubber
- Spent tubes assemblies can be sent to the factory for disposal.

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In addition to the above, mammography tubes assemblies contain an output window made by beryllium, which is highly toxic if ingested.

This part must be in no way machined or worked in a manner suitable of producing powders or fumes, sand blasted or acid etched.

We recommend to send the spent mammography tubes assemblies to the factory, leaving them untouched.

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IAE performs all the possible actions, at design, manufacturing and test level, to make safe the use of the tubes assemblies of its production.

However there exist residual risks connected to the physical principles of operation, which are characterized by high voltages, temperatures and energies being present in the x-ray tube assembly, and by a limited and random life of the components.

These risks are to be evaluated vs. the diagnostic benefits and be well known to those who must take precautions on the subject.

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## 11.1 Radiazioni - Radiations - Radiations

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X-rays emitted by the tube assembly vary in intensity and energy according to the equipment setting.

While being useful from a diagnostic viewpoint, X-rays are anyway harmful to the patient.

The equipment use must thus be strictly limited to qualified personnel under the guide of a radiologist, who can determine the correct exposure parameters and establish a risk to benefit ratio.

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In case of an incorrect assembling of added filtration by the installer, the patient can be exposed to an undue and harmful quantity of low energy radiation.

*The installer must check that the total filtration of the useful beam complies with the existing regulations.*

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In case of incorrect assembling, the tube assembly can give an inadequate protection against leakage radiation.

An excessive leakage radiation exposes the patient to a higher absorbed dose and the operator to a dose in excess of the admissible level.

*The equipment designer must design the leakage radiation shielding so that the shielding continuity between the tube unit and the collimator is insured.*

*The installer must assemble the equipment shielding correctly and check their integrity.*

*The user must have his equipment checked periodically by a qualified expert.*

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The tube assembly is a weight that can be located above the patient.

An accidental drop can cause serious injuries or death of the patient.

*The installer must comply with the assembling instructions of this manual and the equipment manual.*

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The tube assembly and the radiological equipment operate with extremely high voltages.

In normal use the protection from these voltages is ensured by the

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equipment, who are so designed as to prevent any contact with the parts under tension during use.

During installation and maintenance this protection can be removed.

*For this reason the installer and the service engineer must cut off power to the equipment and discharge the cables and the filter capacitors prior than having access to the high voltage terminals.*

*After installation or maintenance, the cable terminals must be tightened so that access to these by the user is made impossible.*

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In case of tube assembly failure or other high voltage parts failure, arcs can occur between the parts connected to high voltage and the tube housing metal envelope that is normally at the ground potential.

The hazard is an electrical shock to the patient, the operator or the installer.

*The equipment designer must design the ground connections in function of the high energies and capacitances involved, and determine the ground paths so that they do not cross the patient and operator contact areas.*

*The installer must take care that the high voltage cables and the ground connections are installed correctly.*

*The operator must suspend the equipment operation as soon as an abnormal operation of high voltage circuits is detected or suspected.*

*The service engineer must check the tube and generator high voltage stability, cables integrity and cable terminals insulation.*

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The tube assembly external surface can reach 80°C.

This temperature is high enough to cause scalding by contact to the patient and the operator.

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The tube assembly contains insulating oil and a glass vacuum tube whose anode reaches high temperatures.

In case of equipment malfunctioning or misuse, the housing can emit high temperature insulating oil, which can cause severe burning to the patient and the operator.

*The designer must comply with the pertinent generators regulations, include connections for all the tube assembly thermal safety devices and design suitable circuits to limit the energy sent to the tube.*

*The installer must connect all the tube assembly safety devices.*

*The user must respect the tube assembly cooling delays.*

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In case of alteration of the focal spot size compared with those required by the pertinent regulations, the image quality can be impaired, causing possible diagnostic errors.

*The user must perform periodical image quality tests, as required by the pertinent regulations.*

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When servicing the equipment for external and internal cleaning of tube assembly, it is possible to come in contact with moving parts or hot surfaces.

*The service engineer must disconnect the equipment from power supply before performing any cleaning.*

*The equipment designer must protect the tube assembly with an adequate casing, so that access from users is prevented.*

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If incorrect power supplies are sent to the device, tube assembly components rupture, arcs towards ground or device operation failure can intervene.

The following incorrect power supplies are evaluated:

Filament overvoltage or overcurrent: filament is melted and stops tube operation without any further damage.

Excessive power to stator: the thermal switch cuts power supply to stator without any damage.

Air circulator overvoltage: air circulator may fail, making air cooling fail and tube thermal switch open.

Anode overvoltage can give rise to a tube arc that is discharged through the ground connection.

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All the seasoning procedure should be performed by the large focal spot. Stator boost cycles must be inhibited during serial exposures in order to prevent tube unit overheating and stator damage.

Perform the following steps by selecting the technique as shown in the tables.

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Perform the required number of exposures with kV, mA, exposure time and cooling delays as shown in the table.

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Starting from the listed initial kV values, make a set of exposures with the required kV, mA, exposure time and cooling delays. At the end of the series increase kV by the required step and repeat the series. When the final kV value is attained, make the final number of exposures as shown in the table.

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Starting from the initial kV value, perform fluoroscopy at the listed kV, mA, fluoroscopy time and cooling delay. Step up kv by the listed value and repeat fluoroscopy. When the final kV value is attained, make the final number of fluoroscopy cycles.

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Starting from the initial kV value, make radiographic exposures at the listed kV, mA, exposure time and cooling delay. At the end of the series increase kV by the required step and repeat the series. When the final kV value is attained, make the final number of exposures as shown in the table.

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If normal operation does not attain the maximum kV values, the seasoning cycles can be limited to the kV values of actual use.  
The warm-up phase is also intended as a daily warm-up.

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The local heat exchanger is supplied as an option to the C52 super housing.

It is a finned copper tube exchanger, forced air cooled by two axial fans. Oil is circulated by a magnetic drive centrifugal pump.

The assembly is contained in a stainless steel casing, and is permanently connected to the housing by two flexible hoses. Maximum hoses length is 1 meter.

The system comes oil filled under vacuum from the factory, and must not be disconnected on field.

*If the installation requires longer oil hoses, or the capability of disconnecting and reconnecting hoses, a remote heat exchanger, not described in this section, can be supplied.*

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During any system handling, from unpacking up to tube installation, never grip the parts by the oil hoses, and avoid any traction, torque or bending stress on hoses.

Since the system is sealed under vacuum, oil hoses and couplings must not be disconnected.

Check the oil system for shipping damage. If any part appears to be loose, cracked or leaks oil, air may have entered the system. Suspend installation and contact the factory if this is the case.

Install the housing according to instructions in par. 7. Heat exchanger casing should be fastened by adequate mechanical means, complying to regulations concerning suspended masses under motion.

The casing contains electrical parts under mains voltage, and must be protected against liquid spills.

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An air free path having the same cross section as the fans window must be ensured.

Connect the heat exchanger to power supply at the nominal voltage by a cable type IEC/CEE22 with ground connection. Comply with local regulations pertaining to wires gauge, insulation and grounding.

Safety against possible cooling system malfunctioning is ensured by tube housing thermal switches.

To protect patients and personnel against personal injury, it is mandatory to connect these safety devices to the equipment, and verify that an adequate reaction from the equipment is obtained by their intervention.

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No user serviceable parts are contained in the system.

Power supply socket contains a fuse that must be replaced, if needed, by a fuse with same rating, after having found and removed the cause of failure.

Heat exchanger fins must be kept free from dust in order to maintain an adequate air flow. Accumulated dust should be removed by a vacuum cleaner and air blast. The grids on the fins side only may be removed for this purpose, but care must be taken in avoiding damages to the fins.

Grids on the fans side must not be removed, because inner parts would fall inside.

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