


|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 19 |
|  |  | Page    | 1/137        |

# Mars1717V User's Manual



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 3/137        |

## To Customers



Congratulations on your purchase of the Mars1717V Wireless Digital Flat Panel

(hereinafter referred to as Mars1717V) which is manufactured by iRay Technology (Shanghai) Ltd. (hereinafter referred to as iRay).

At iRay, we strive to not only make the world-class products that deliver the best value possible to our customers but also offer the highest quality of service and customer care. Please take time to read through this user guide in order to utilize the product effectively. We hope you enjoy the experience with iRay Mars1717V.

If you have any questions or suggestions, please feel free to contact us.

### Service Office

**Tel: +86 21 50720560 - 8059**

**Fax: +86 21 50720561**

**E-mail: service@iraychina.com**

**Location: 2F, Building 7, No.590, Ruiqing Rd, Pudong, Shanghai, China PC: 201201**

## Notes on usage and management of the equipment

- Read all of the instructions in the user guide before your operation. Give particular attention to all safety precautions.
- Only a physician or a legally certified operator should use this product.
- The equipment should be maintained in a safe and operable condition by maintenance personnel.
- Use only computers and image display monitors complying with IEC 60601-1 or IEC 60950-1. For details, consult our sales representative or local iRay dealer.
- Use only the dedicated cables. Do not use any cables other than those supplied with this product.
- Request your sales representative or local iRay dealer to install this product

## Caring for your environment



This symbol indicates that this product is not to be disposed of with your residential or commercial waste.

## Recycling iRay Equipment

Please do not dispose of this product with your residential or commercial waste. Improper handling of this type of waste could have a negative impact on health and on the environment. Some countries or regions, such as the European Union, have set up systems to collect and recycle electrical or electronic waste items.

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 4/137        |

Contact your local authorities for information about practices established in your region. If collection systems are not available, call iRay Customer Service for assistance.

## Disclaimer

1. iRay shall not be liable to the purchaser of this product or third parties for any damage, losse, or injury incurred by purchaser or third parties as a result of fire, earthquake, any accident, misuse or abuse of this product.
2. iRay shall not be liable to any damage, loss, or injury arising from unauthorized modifications, repairs, or alterations to this product or failure to strictly comply with iRay's operating and maintenance instructions.
3. iRay shall not be liable for any damage or loss arising from the use of any options or consumable products other than those dedicated as Original iRay Products by iRay Technology.
4. It is the responsibilities of the user/attending physicians for maintaining the privacy of image data and providing medical care services. iRay shall not be responsible for the legality of image processing , reading and storage nor it shall be responsible for loss of image data for any reason.
5. Information regarding specification, compositions, and appearance of this product is subject to change without prior notice.

## Copyright

- All rights reserved
- No part of this publication may be reproduced in any form or by any means without the written permission of iRay. The information contained herein is designed only for use with iRay Mars1717V.








## Trademarks

- The iRay name and iRay logo are registered trademarks of iRay Technology (Shanghai) Ltd.

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 5/137        |











## Symbols and Conventions

The following symbols and conventions are used throughout the user guide.









|   |   |
|---|---|
|    | <p>This symbol is used to identify conditions under which improper use of the product may cause death or serious personal injury.</p> |
|    | <p>This notice is used to identify conditions under which improper use of the product may cause minor personal injury.</p>            |
|    | <p>This notice is used to identify conditions under which improper use of the product may cause property damage.</p>                  |
|  | <p>This is used to indicate a prohibited operation.</p>   |
|  | <p>This is used to indicate an action that must be performed.</p>   |
|  | <p>This is used to indicate important operations and restrictions.</p>  |
|  | <p>This is used to indicate operations for reference and complementary information.</p>   |



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 7/137        |

|   |  |
|---|--|
|    | This symbol is used to indicate the name and address of the manufacturer.                                      |
|    | This symbol is used to indicate the name and address of iRay authorized representative in the European region. |
|    | This symbol is used to indicate consultation of the user guide for general information.                        |
|    | Safety Signs: please refer to the user guide for safety instructions   |
|   | Safety Signs: Dangerous Voltage  |
|  | Stand-by   |
|  | Handled with care  |
|  | FPD is allowed to withstand 100 kg on the surface  |
|  | This symbol is used to indicate the operational temperature limits.  |
|  | This symbol is used to indicate the storage temperature limits.  |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 8/137        |

|   |  |
|---|--|
|    | non-ionizing radiation   |
| <b>FCC</b>  | Federal Communications Commission certificate                          |
|    | Package symbol, fragile.   |
|    | Package symbol, keep away from sunlight                                |
|    | Package symbol, keep dry   |
|  | Package symbol, this symbol is used to indicate the humidity limits.   |
|  | Package symbol, keep the equipment up right                            |
|  | Package symbol, do not roll the transportation package.                |
|  | Package symbol, this symbol is used to indicate stacking limit number. |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 9/137        |

# Contents

|   |           |
|---|-----------|
| <b>To Customers</b>                         | <b>3</b>  |
| <b>Symbols and Conventions</b>              | <b>5</b>  |
| <b>Labels and markings on the equipment</b> | <b>6</b>  |
| <b>Contents</b>                             | <b>9</b>  |
| <b>1 Safety Information</b>                 | <b>13</b> |
| 1.1 Safety precautions                      | 13        |
| 1.2 Notes for Using                         | 18        |
| <b>2 General Description</b>                | <b>20</b> |
| 2.1 Scope                                   | 20        |
| 2.2 Model                                   | 21        |
| 2.3 Characteristic                          | 21        |
| 2.4 Intended use                            | 21        |
| 2.5 Standard Product Components             | 22        |
| 2.6 Optional Product Component              | 24        |
| 2.7 Components Description                  | 25        |
| 2.8 Product Specification                   | 28        |
| 2.8.1 Detector                              | 28        |
| 2.8.2 Battery                               | 30        |
| 2.8.3 Battery Charger                       | 31        |
| 2.8.4 Power supply                          | 32        |
| 2.8.5 Infrared Device (Optional)            | 32        |
| 2.8.6 AP Router (Optional)                  | 32        |
| 2.8.7 Wireless Communication                | 32        |
| 2.8.8 Recommended Appliance Condition       | 33        |
| 2.8.9 Mechanical Outlines                   | 34        |
| 2.8.10 Use Environment                      | 36        |
| <b>3 Installation</b>                       | <b>37</b> |
| 3.1 Detector Installation                   | 37        |
| 3.2 Battery Charging Installation           | 40        |



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 10/137       |

|          |   |           |
|----------|---|-----------|
| 3.3      | Software Installation                               | 41        |
| 3.4      | Panel Infrastructure                                | 41        |
| <b>4</b> | <b>Operation</b>                                    | <b>62</b> |
| 4.1      | Main Operation                                      | 62        |
| 4.1.1    | Software Mode                                       | 62        |
| 4.1.2    | Inner Mode  | 65        |
| 4.1.3    | Isync Plus Mode                                     | 67        |
| 4.2      | Connection Build                                    | 69        |
| 4.3      | Panel Configuration                                 | 70        |
| 4.4      | Correction and Calibration Template Generation      | 70        |
| 4.4.1    | Pre-offset Template Generation                      | 70        |
| 4.4.2    | Gain Calibration Template Generation                | 71        |
| 4.4.3    | Defect Correction Template Generation               | 73        |
| 4.4.4    | Lag Template Generation                             | 76        |
| 4.5      | Image Acquisition Continually                       | 77        |
| 4.6      | Image save  | 77        |
| 4.6.1    | Multiple images                                     | 77        |
| 4.6.2    | Single image  | 78        |
| 4.7      | Image Check and upload                              | 79        |
| 4.7.1    | Local Image Check                                   | 79        |
| 4.7.2    | Panel Image Upload                                  | 80        |
| 4.7.3    | Panel Image Check                                   | 83        |
| 4.8      | Defect Template Check and Modification              | 83        |
| 4.8.1    | Defect Template Check                               | 83        |
| 4.8.2    | Defect Template Modification                        | 84        |
| 4.9      | Correction and Calibration Management               | 86        |
| 4.9.1    | Correction and Calibration template synchronization | 86        |
| 4.9.2    | Correction and Calibration management               | 88        |
| 4.10     | Sleep and Wake Up                                   | 90        |
| 4.10.1   | Sleep   | 91        |
| 4.10.2   | Wake Up   | 93        |
| 4.11     | Firmware Update                                     | 93        |
| 4.11.1   | MCU Update  | 93        |
| 4.11.2   | FTP Server  | 95        |
| 4.11.3   | Firmware Update                                     | 97        |
| 4.12     | Short cut   | 99        |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 11/137       |

|          |  |            |
|----------|--|------------|
| 4.13     | Software                                       | 100        |
| 4.13.1   | Main GUI                                       | 100        |
| 4.13.2   | Message Box                                    | 101        |
| 4.13.3   | Configuration GUI                              | 103        |
| 4.13.4   | Correction and Calibration                     | 110        |
| 4.13.5   | Image Check and upload                         | 112        |
| 4.13.6   | Infrared Registration                          | 113        |
| <b>5</b> | <b>Regulatory Information</b>                  | <b>115</b> |
| 5.1      | Medical equipment safety standards             | 115        |
| 5.2      | Guidance and manufacture's declaration for EMC | 116        |
| 5.3      | Radio Frequency Compliance Information         | 120        |
| 5.3.1    | FCC Compliance                                 | 121        |
| 5.4      | Battery Safety Standards                       | 121        |
| 5.5      | Product Label                                  | 122        |
| 5.5.1    | Detector                                       | 122        |
| 5.5.2    | Battery  | 124        |
| 5.5.3    | Battery Charger                                | 125        |
| <b>6</b> | <b>Trouble Shooting</b>                        | <b>126</b> |
| 6.1      | Cable Inspection                               | 126        |
| 6.1.1    | Ethernet cable                                 | 126        |
| 6.1.2    | Date cable inspection                          | 127        |
| 6.2      | Adapter Inspection                             | 127        |
| 6.2.1    | Adapter connector definition                   | 127        |
| 6.2.2    | Inspection method                              | 128        |
| 6.3      | Battery Inspection                             | 128        |
| 6.3.1    | Battery pin definition                         | 128        |
| 6.3.2    | Inspect Method                                 | 129        |
| 6.4      | Dock Inspection                                | 129        |
| 6.4.1    | Dock pin definition                            | 129        |
| 6.4.2    | Inspection method                              | 130        |
| 6.5      | Detector Main Unit Inspection                  | 130        |
| 6.5.1    | Get SDK and detector log                       | 130        |
| 6.6      | Power up failure inspection                    | 132        |
| 6.7      | Connection inspection                          | 132        |
| 6.7.1    | Connection failed in wired mode                | 132        |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 12/137       |




|          |   |            |
|----------|---|------------|
| 6.7.2    | Connection Failed in Client Mode                    | 133        |
| 6.7.3    | Disconnect in AP Mode                               | 133        |
| 6.8      | Dark Image Acquisition After Exposure in Inner Mode | 133        |
| 6.9      | No Image Acquire after Exposure                     | 134        |
| 6.10     | Image Inspection                                    | 134        |
| 6.10.1   | Artifact on bright image                            | 134        |
| 6.10.2   | Horizontal Line on bright image                     | 134        |
| 6.10.3   | Defective line                                      | 136        |
| 6.10.4   | Defective point                                     | 136        |
| 6.10.5   | TFT broken  | 136        |
| <b>7</b> | <b>Service Information</b>                          | <b>137</b> |
| 7.1      | Product lifetime                                    | 137        |
| 7.2      | Regular inspection and Maintenance                  | 137        |
| 7.3      | Repair  | 138        |
| 7.4      | Replacement parts support                           | 138        |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 13/137       |

# 1 Safety Information

## 1.1 Safety precautions

Follow these safeguards and properly use the equipment to prevent injury and damage to any equipment/data.

| WARNING   |  |
|---|--|
| <b>Installation and environment of use</b>  |  |
|    | <ul style="list-style-type: none"> <li>Do not use or store the equipment near flammable chemicals such as alcohol, thinner, benzene, etc.</li> </ul> <p>If chemicals are spilled or evaporate, it may result in fire or electric shock through contact with electric parts inside the equipment. Also, some disinfectants are flammable. Be sure to take care when using them.</p> |
|    | <ul style="list-style-type: none"> <li>Do not connect the equipment with anything other than specified.</li> </ul> <p>Doing so may result in fire or electric shock.</p>   |
|   | <ul style="list-style-type: none"> <li>All the patients with active implantable medical devices should be kept away from the equipment.</li> </ul>   |
| <b>Power supply</b>   |  |
|  | <ul style="list-style-type: none"> <li>Do not operate the equipment using any type of power supply other than the one indicated on the rating label.</li> </ul> <p>Otherwise, it may result in fire or electric shock.</p>   |
|   | <ul style="list-style-type: none"> <li>Do not handle the equipment with wet hands.</li> </ul> <p>You may experience electric shock that could result in death or serious injury.</p>   |
|   | <ul style="list-style-type: none"> <li>Do not place heavy object such as medical equipment on cables and cords. Do not pull, bend, bundle, or step on them to prevent their sheath from being damaged, and do not alter them neither.</li> </ul> <p>Doing so may damage the cords which could result in fire or electric shock.</p>  |
|   | <ul style="list-style-type: none"> <li>Do not supply power to more than one piece of equipment using the same AC outlet.</li> </ul> <p>Doing so may result in fire or electric shock.</p>  |
|   | <ul style="list-style-type: none"> <li>Do not turn ON the system power when condensation has formed on the equipment.</li> </ul> <p>Doing so may result in fire or electric shock.</p>   |
|   | <ul style="list-style-type: none"> <li>Do not connect a multiple portable socket-outlet or extension cord to the system.</li> </ul>  |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 14/137       |

Doing so may result in fire or electric shock.

- **To avoid the risk of electric shock, this equipment must only be connected to power supply with protective earth.**

Not doing so may result in fire or electric shock.



- **Securely plug the power cord into the AC outlet.**

If contact failure occurs, or if metal objects come into contact with the exposed metal prongs of the plug, fire or electric shock may result.

- **Be sure to turn OFF the power to each piece of equipment before connecting or disconnecting the cords.**

Otherwise, you may get an electric shock that could result in death or serious injury.

- **Be sure to hold the plug or connector to disconnect the cord.**

If you pull the cord, the core wire may be damaged, resulting in fire or electric shock.

## WARNING

### Handling



- **Never disassemble or modify the equipment. No modification of this equipment is allowed. Parts of the Mars1717V that are not serviced or maintained while in use with the patient**

Doing so may result in fire or electric shock. Also, since the equipment incorporates parts that may cause electric shock as well as other hazardous parts, touching them may cause death or serious injury.

- **Do not place anything on top of the equipment.**

The object may fall and cause an injury. Also, if metal objects such as needles or clips fall into the equipment, or if liquid is spilled, it may result in fire or electric shock.

- **Do not hit or drop the equipment.**

The equipment may be damaged if it receives a strong jolt, which may result in fire or electric shock if the equipment is used without being repaired.

- **Do not put the equipment and pointed objects together.**

The equipment may be damaged. If so, the equipment should be used in bucky.



- **Have the patient take a fixed posture and do not let the patient touch parts unnecessarily.**

If the patient touches connectors or switches, it may result in electric shock or malfunction of the equipment.

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 15/137       |

### When a problem occurs



- **Should any of the following occurs, immediately unplug the power cord of adaptor or battery, and contact your sales representative or local iRay dealer:**

When there is smoke, an odd smell or abnormal sound.

When liquid has been spilled into the equipment or a metal object has entered through an opening.

When the equipment has been dropped and damaged.

### Maintenance and inspection



Prohibited

- **Please turn OFF the power of the equipment and unplug the power cord of adaptor before cleaning.**
- **NEVER use alcohol, ether and other flammable cleaning agent for safety. NEVER use methanol, benzene, acid and base because they will erode the equipment.**
- **DON'T dip the equipment into the liquid.**
- **Please make sure that the equipment's surface & plugs are dry before turning ON.**

Otherwise, it may result in fire or electric shock.



- **Clean the plug of the power cord periodically by unplugging it from the AC outlet and removing dust or dirt from the plug, its periphery and AC outlet with a dry cloth.**

If the cord is kept plugged in for a long time in a dusty, humid or sooty place, dust around the plug will attract moisture; this could cause insulation failure that may result in a fire.

- **For safety reasons, be sure to turn OFF the power to each piece of equipment when performing inspections indicated in this manual.**

Otherwise, electric shocks may occur.

## CAUTION

### Installation and environment of use



- **Do not install the equipment in any of the locations listed below. Doing so may result in failure, malfunction, equipment falling, fire or injury.**

Close to facilities where water is used

Where it will be exposed to direct sunlight

Close to the air outlet of an air-conditioner or ventilation equipment

Close to heat source such as a heater

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 16/137       |

Where the power supply is unstable

In a dusty environment

In a saline or sulfurous environment

Where temperature or humidity is high

Where there is freezing or condensation


In areas prone to vibration

On an incline or in an unstable area

- **Take care that cables do not become tangled during use. Also, be careful not to get your feet caught by cable.**

Otherwise, it may cause a malfunction of the equipment or injury of the user due to tripping over the cable.

- **Non-medical equipment such as the battery charger, access point and IR data communication unit cannot be used in patient's vicinity.**



## Power supply



- **Always connect the three-core power cord plug to a grounded AC power outlet.**
- **To make it easy to disconnect the plug at any time, avoid putting any obstacles near the outlet. Otherwise, it may not be possible to disconnect the plug in an emergency.**
- **Be sure to ground the equipment to an indoor grounded connector. Also, be sure to connect all the grounds for the system to a common ground.**
- **Do not use any power source other than the one provided with this equipment.**

Otherwise, fire or electric shock may be caused due to leakage.

## Handling



- **Do not spill liquid or chemicals onto the equipment. In case the patient is injured, it is not allowed to contact with blood or other body fluids.**

Doing so may result in fire or electric shock.

In such a situation, protect the equipment with a disposable cover as necessary.

- **Turn OFF the power and pull out the plug to each piece of equipment for safety when**

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 17/137       |


not used.

## CAUTION

### Handling



- **Handle the equipment carefully.**
- **Do not submerge the equipment in water.**
- **The internal image sensor may be damaged if something hits against it or it is dropped.**  
If the equipment is dropped, the drop sensor inside will turn red and the equipment will not be warranted by iRay.





- **Do not place excessive weight on the equipment.**  
Otherwise, the internal image sensor may be damaged and image may be incorrect.

<Load Limit>

Uniform load: 100 kg over the whole area of the detector surface.


Local load: 100 kg on an area 4 cm diameter.



- **Be sure to use the equipment on a flat surface so it will not bend. Otherwise, the internal image sensor may be damaged. Be sure to securely hold the detector while using it in upright positions. Otherwise, the detector may fall over, resulting in injury to the user or patient, or may flip over, resulting in damage to the inner device.**



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 18/137       |



Keep the same load (same pressure) on the detector when acquiring the image. Or the image will be incorrect.

## CAUTION



- Do not close to fire, do not use in high temperature
- Do not invert positive and negative pole
- Do not contact with metal in case of short circuit
- Do not insert sharp objects into battery
- Do not beat battery
- Do not stand on battery
- Do not use battery out of rules
- Do not dispose battery or change internal structure
- Do not submerge battery in water, please keep dry in storage and do not contact with water in use
- Please charge battery with charger following IEC60601-1 & IEC62133 Standards provide by us
- Do not mix battery with ones not provided by our company
- Do not charge battery with broken charger.

## 1.2 Notes for Using

When using the equipment, take the following precautions. Otherwise, problems may occur and the equipment may not function correctly.

### Before exposure

- Be sure to check the equipment daily and confirm that it works properly.
- Be sure there be a battery installing on the Mars1717V to avoid the power off suddenly.
- Sudden heating of the room in cold areas will cause condensation to form on the equipment. In

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 19/137       |

this case, wait until the condensation evaporates before performing an exposure. If the equipment is used while condensation is formed on it, problems may occur in the quality of captured images. When an air-conditioner is used, be sure to raise/lower the temperature gradually so that a difference of temperature in the room and equipment does not occur, to prevent condensation.

- The detector should warm up for 15 minutes before exposure or updating the gain map or defect map.

#### **During exposure**

- Do not move the power or Ethernet Cables during exposure, or it may cause image noise or artifacts, even incorrect images.
- Do not use the devices near the equipment generating a strong magnetic field. Otherwise, it may cause image noise, artifacts or even incorrect images.

#### **Disinfection and Cleaning**

- After every examination, wipe the patient contact surfaces of the detector using disinfectants such as ethanol, to prevent the risk of infection. For details on how to sterilize, consult a specialist.
- Do not spray the detector directly with disinfectants or detergents.
- Wipe it with a cloth slightly damped with a neutral detergent. Do not use solvents such as alcohol, thinner, benzene, acid and base. Doing so may damage the surface of the equipment.
- It's recommended to use a waterproof non-woven cover as the isolated layer between detector and the bleeding patient.


|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 20/137       |

## 2 General Description

Mars1717V is a cassette-size wireless X-ray flat panel detector based on amorphous silicon thin-film transistor technologies. It is developed to provide the highest quality of radiographic image, which contains an active matrix of 3072×3072 with 139um pixel pitch. Detectors' scintinator has two options which are Standard GOS(Gadolinium Sulfoxylate) and CsI(CaesiumIodide). However the most great improvement is Mars1717V supports wireless communication between detectors and PC. Mars1717V's power supply includes battery. Mars1717V can be used as a real portable panel.


### 2.1 Scope

This manual contains information about the iRay Mars1717V. All operators must read and understand this manual before using equipment. All information in this manual, including the illustrations, is based on equipment prototype. If configuration of your equipment does not have any of these items, information about these items in the manual does not apply to your equipment.



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 21/137       |

## 2.2 Model



Product Type: Battery-KV-----Rechargeable lithium battery

Product Type: Charger-KV-----Battery charger

## 2.3 Characteristic

- Wireless static Flat Panel Detector used for general radiography.
- 17 × 17 inch
- Sync-Shot exposure trigger
- GOS or CsI scintillation screen.
- Easy to change the cable and upgrade firmware.
- Battery recycling

## 2.4 Intended use

This equipment provides digital X-ray imaging for diagnosis of disease, injury, or any applicable health problem. The image is obtained as the result of X-ray passing through the human body and detected by the equipment. This device is intended to be used in the holder or bucky which is well insulated to panel. The holder or the bucky is well grounded.

iRay will provide equipment and software support for integration of system. The length of both Ethernet Cable and DC Power Cable cannot exceed 3.5 m. or the impedance of protective earth connections may exceed the safety threshold.

This panel is not intended for mammography or dental applications.

According to the Mars1717V intended use and the result of risk management, identifying and describing the essential performance as the following:


- To get imaging of dark field, the Mars1717V shall be not influenced to the imaging acquisition

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 22/137       |

- b) To keep the data transmission function, the Mars1717V shall be not influenced to the data and signal transmission



## 2.5 The relative position between patient and detector

Because of the crosstalk effect of Amorphous silicon flat-panel detector, Pay attention to the relative position of patient and detector, the recommended position as shown below, Otherwise, the image is prone to abnormal light lines.







## 2.6 Standard Product Components

Mars1717V comes with power supply both 24V DC and battery package. Once powered on, it would build a connection with PC through Ethernet cable or Wireless connection.

|  | Item   | Description        |
|--|--|--------------------|
| Mars1717V Detector   |  | 1 pcs<br>Main Unit |
| Medical Adapter for<br>• Detector and<br>• Battery Charger |  | 2 pcs<br>DC 24V    |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 23/137       |

|                        |  |                     |
|------------------------|--|---------------------|
| Battery Pack           |    | 2 pcs<br>Battery-KV |
| Ethernet Cable         |    | 1 pcs<br>3.5 m      |
| Gigabit Ethernet Cable |   | 1 pcs<br>3 m        |
| AC Power Cable         |  | 2 pcs               |
| DC Power cable         |   | 1 pcs<br>3.5 m      |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 24/137       |

|                 |  |   |
|-----------------|--|---|
| Battery Charger |  | 1 pcs   |
| CD-Rom          |   | 1 pcs<br>Gain correction data<br>Defect correction map<br>SDK<br>Manual |


## 2.7 Optional Product Component

|                    | Item   | Description |
|--------------------|--|-------------|
| Wireless AP Device |   | 1 pcs       |
| Infrared Device    |  | 1 pcs       |


|  |                                 |         |            |
|--|---------------------------------|---------|------------|
|  | iRay Technology (Shanghai) Ltd. | Doc NO. | 037-201-02 |
|  |                                 | Version | A2         |
|  | Mars1717V User' s Manual        | Date    | 2016.02.23 |
|  |                                 | Page    | 25/137     |

## 2.8 Components Description


### 2.8.1.1 Detector



#### External Signals Input



#### Control Panel



#### Control Panel


| Item | Name               | Description                         |
|------|--------------------|-------------------------------------|
| A    | DC Input Interface | 24V DC input                        |
| B    | Ethernet Interface | Gigabit Ethernet Wire               |
| C    | Detector Indicator | Detector indicator of control panel |
| D    | Power Button       | Power button of control panel       |



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 26/137       |

|   |                   |  |
|---|-------------------|--|
| E | Infrared Window   | Infrared device window                   |
| F | Antenna           | Antenna                                  |
| G | Handle            | For the panel carried                    |
| H | Maintenance Cover | For service engineer to maintenance      |
| I | Battery Lock      | The lock button for detaching battery    |
| J | Detector Label    | Product information.                     |
| K | Switch            | AP/Client Mode Switch                    |
| L | Carbon Film       | Panel Carbon Film, have biocompatibility |

### 2.8.1.2 Battery



| Item | Name              | Description                      |
|------|-------------------|----------------------------------|
| A    | Battery Label     | /                                |
| B    | Battery Interface | 8 Pin Battery connector          |
| C    | The location pin  | /                                |
| D    | Indicator         | Installation direction indicator |



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 27/137       |

### 2.8.1.3 Battery Charger




| Item | Name                | Description                           |
|------|---------------------|---------------------------------------|
| A    | Battery Interface   | 8 Pin Battery connector               |
| B    | Capacity Indicator  | The indicator definition is as follow |
| C    | Power Indicator     | The indicator definition is as follow |
| D    | Hand Pull Position  | /                                     |
| E    | The limit ball plug | /                                     |
| F    | DC Jack             | 24V DC input                          |





Power indicator definition:

| Power Indicator | Lighting Status   | Operating Status             |
|-----------------|---|------------------------------|
| OFF             |  | No external DC adaptor input |
| GREEN ON        |  | External DC adaptor input    |

The battery charging capacity indicator definition:

| X Group Indicator      | Lighting Status   | Operating Status  |
|------------------------|---|-------------------|
| I, II and III grid OFF |  | No battery Insert |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 28/137       |

|  |   |  |
|--|---|--|
| I grid blinking<br>II and III grid OFF |  | Battery Insert with capacity $\leq 30\%$ , charging  |
| II grid blinking<br>I and III grid OFF |  | Battery Insert with capacity $>30\%$ and $\leq 60\%$ , charging  |
| III grid blinking<br>I and II grid OFF |  | Battery Insert with capacity $>60\%$ and $\leq 95\%$ , charging  |
| I and II grid OFF<br>III grid ON       |  | Battery Insert with capacity $>95\%$ and charging, when capacity = 100%, charging can stop automatically |

#### 2.8.1.4 Power Supply

Mars1717V supports both DC Power and Battery package input.

#### 2.8.1.5 Infrared Device

Mars1717V does not include Infrared Device. Users can select by themselves, however some basic requirements should be followed.

## 2.9 Product Specification

### 2.9.1 Detector

#### 2.9.1.1 Basic

| Item                   | Specification                              |
|------------------------|--|
| Model                  | Mars1717V-PSI (GOS)<br>Mars1717V-VSI (CsI) |
| Image Sensor           | a-Si (Amorphous Silicon) TFT               |
| Pixel Size             | 139 $\mu\text{m}$                          |
| Effective Array        | 3072 x 3072                                |
| Effective Area (H x V) | 427 x 427 mm                               |
| Fill Factor            | 60%  |
| Greyscales             | 14bit                                      |
| Spatial Resolution     | 2.8 Lp/mm (Standard GOS)                   |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 29/137       |

|                                   |  |
|-----------------------------------|--|
|                                   | 3.1Lp/mm (CsI)   |
| Image Acquisition Time (Wired)    | Preview Acquisition Time : 6 sec<br>Processed Acquisition Time : 10 sec. (including Preview Time)                        |
| Image Acquisition Time (Wireless) | Preview Acquisition Time : 9 sec.  |
| Both AP mode and Client mode      | Processed Acquisition Time : 16 sec. (including Preview Time)  |
| Cycle Time                        | Min. 13s @Wired;Min.16s @Wireless  |
| Power Consumption                 | Max. 16W   |
| Dimension (L × W × H)             | 460 x 460 x 15.2 mm  |
| Weight (with one battery)         | Mars1717V-PSI: 4.3 kg without battery, 4.6 kg with battery<br>Mars1717V-VSI: 4.5 kg without battery, 4.8 kg with battery |
| Image Transfer                    | Wired : Gigabit Ethernet(1000BASE-T)<br>Wireless : IEEE802.11a/b/g/n   |
| Data Transmission Rate (Wireless) | 802.11b : Max. 11Mbps<br>802.11a/g : Max. 54Mbps<br>802.11n : Max. 300Mbps (MIMO 2x2)                                    |
| X-ray Energy                      | 40kV to 150kV  |

#### 2.9.1.2 MTF

The MTF with GOS should meet the following table

| Spatial frequency (lp/mm) | MTF (GOS) | MTF(CsI) |
|---------------------------|-----------|----------|
| 1.0                       | 0.39      | 0.60     |
| 2.0                       | 0.12      | 0.31     |
| 3.0                       | 0.04      | 0.15     |


#### 2.9.1.3 DQE

The DQE with GOS should meet the following table

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 30/137       |

| Spatial frequency (lp/mm)@RQA5 | DQE (GOS)@3.2uGy | DQE (CsI)@2.5uGy |
|--------------------------------|------------------|------------------|
| 0                              | 0.30             | 0.52             |
| 1.0                            | 0.15             | 0.39             |
| 2.0                            | 0.05             | 0.24             |
| 3.0                            | 0.01             | 0.16             |

## 2.9.2 Battery



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 31/137       |

| Item                   | Specifications                              |
|------------------------|---|
| Model                  | Battery-KV                                  |
| Rated Capacity         | Min. 3950mAh, Typ. 4180mAh @ Discharge 0.2C |
| Nominal Voltage        | 10.8V                                       |
| Charge Voltage         | 12.6±0.05V                                  |
| Discharged End Voltage | 8.25V                                       |
| Charging Method        | CC-CV                                       |
| Operating Temperature  | Charge 0℃-+45℃, Discharge -20℃-+60℃         |
| Storage Temperature    | 1 month -20℃-+45℃                           |
|                        | 3 month 0℃-+30℃                             |
|                        | 6 month 5℃-+20℃                             |
| Relative Humidity      | 65±20%                                      |
| Dimension (L × W × H)  | 210 x 115 x 7.5 mm                          |
| Weight                 | 0.29 kg                                     |

### 2.9.3 Battery Charger



| Item                  | Specifications    |
|-----------------------|-------------------|
| Model                 | Charger-KV        |
| Simultaneous Charging | 2 battery packs   |
| Full charging time    | 2 hours           |
| Rated power supply    | 24V(DC)           |
| Dimension (L × W × H) | 300 x 263 x 42 mm |
| Weight                | 1.26 kg           |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 32/137       |

## 2.9.4 Power supply

Mars1717V supports both DC Power and Battery package input.

| Item            | Specifications       |
|-----------------|----------------------|
| <b>DC Power</b> | <b>24V(DC), 0.8A</b> |
| Battery Package | 10.8V(DC),1.78A      |

## 2.9.5 Infrared Device (Optional)

Mars1717V does not include Infrared Device. Users can select by themselves, however some basic requirements should be followed.

| Item          | Specifications                    |
|---------------|-----------------------------------|
| IRDA Protocol | Compliant with IrDA V1.0 and V1.1 |
| USB           | Compliant with USB V2.0 and V1.1  |
| Data Rate     | Max. 4Mbps                        |

## 2.9.6 AP Router (Optional)

Mars1717V do not include AP Router. Users can select AP Router as they wish, however specification below is a requirements.

| Item               | Specifications  |
|--------------------|---|
| Wireless Standard  | IEEE 802.11 a/b/g/n   |
| Frequency Range    | 2.412 ~ 2.4835 GHz and 5.15 ~ 5.85 GHz  |
| Wireless Data Rate | 802.11b : Max. 11Mbps<br>802.11a/g : Max. 54Mbps<br>802.11n : Max. 300Mbps (MIMO 2x2) |
| Wired Data Rate    | Max. 1Gbps  |

## 2.9.7 Wireless Communication

| Item              | Description   |
|-------------------|---|
| Wireless Standard | IEEE802.1a/b/g/n  |
| Frequency Range   | 2.4G: 2.412 ~ 2.4835 GHz<br>11: (Ch. 1-11) – United States<br>13: (Ch. 1-13) – Europe<br>14: (Ch. 1-14) – Japan |


|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 33/137       |

|                        |   |
|------------------------|---|
|                        | 5G: 5.15 ~ 5.85 GHz<br>12: United States<br>19: Europe<br>8: Japan  |
| Data Transmission Rate | 802.11b : Max. 11Mbps<br>802.11a/g : Max. 54Mbps<br>802.11n : Max. 300Mbps (MIMO 2x2)                               |
| Modulation             | 802.11b:<br>CCK, DQPSK, DBPSK<br>802.11a/g:<br>64 QAM, 16 QAM, QPSK, BPSK<br>802.11n:<br>64 QAM, 16 QAM, QPSK, BPSK |
| Transmission Power     | Max.17dBm   |
| Security               | WPA, WPA-PSK, WPA2, WPA2-PSK, WEP 64bit & 128bit  |
| Antenna                | 2 Dual Band internal antenna  |


### 2.9.8 Recommended Appliance Condition


| Item             | Description  |
|------------------|--|
| Operating System | Windows XP/7 32/64bit  |
| CPU              | Intel Core i73.6G  |
| Memory           | 4G DDR3  |
| Hard Disk        | 160 G  |
| LAN Card         | Intel Pro EXP9301CT PRO<br>Gigabit Network Adapter with PCIe interface |







|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 34/137       |

## 2.9.9 Mechanical Outlines



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 35/137       |



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 36/137       |

## 2.9.10 Use Environment

|  | Temperature | Temperature change | Humidity  | Atmospheric Pressure    | Pressure Change                |
|--|-------------|--------------------|-----------|-------------------------|--------------------------------|
| Operating  | 5~30℃       | <1k/min            | 45~85% RH | 700~1000hPa<br>0~70% RH | <10kp/min<br>(1kp=1.0197E-5Pa) |
| Storage  | -20~60℃     | <1k/min            | 45~85%RH  | 700~1000hPa<br>0~70% RH | <10kp/min<br>(1kp=1.0197E-5Pa) |
| The Mars1717V serial detectors shall operate at an altitude specified not more than 3000m. |             |                    |           |                         |                                |




|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 37/137       |

## 3 Installation

### 3.1 Detector Installation

#### 3.1.1.1 Attach Battery Pack




Mars1717V can be powered on by both battery package and DC power. Once battery package is inserted or DC power is on, Detectors would be activated immediately. If none of battery and DC power is on, Mars1717V would power off. Please see below for battery installation.

|  |  |
|--|--|
| Make sure that the connectors on the battery package are pointed to the cave in battery compartment. |   |
| Slide battery package into battery compartment<br>(Make sure battery capacity overpass 10%) .        |  |
| Slide the battery lock lever.  |  |

#### 3.1.1.2 Attach DC Power


Please see below for DC power installation.

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 38/137       |

|  |   |
|--|---|
| <p>Connect one end of DC Power Cable to the Medical Adapter</p>  |   |
| <p>As figure is power interface and Ethernet interface</p>   |   |
| <p>Connect another end of DC Power Cable to the DC input of the detector.</p> <p>If user want to use wired Ethernet, connect the Ethernet cable to Ethernet interface.</p> |  |

### 3.1.1.3 Booting Up

On the control panel, user can press power button to power on/off.




|   |  |
|---|--|
| <p>If detector is powered off, users can press the button for 4 seconds to power on the detector when battery is inserted and battery capacitor is not less than 10%, or direct current power is connected.</p> <p>If detector is powered on, users can press the button for 4 seconds to shut down the detector. On the other hand, it can also be used as reset internal control IC when button is active for 8s.</p> |  |
|---|--|

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 39/137       |



After booting up the detector, use can check the status LED indicator of detector.

| Power Indicator                 | Lighting Status   | Operating Status |   |          |
|---------------------------------|---|------------------|---|----------|
|                                 |   | Operating        | Battery Capacity  | DC Input |
| OFF                             |    | Power OFF        | /   | /        |
| Orange ON                       |    | Power ON         | ≤10%  | NO       |
| Green ON                        |    | Power ON         | <ul style="list-style-type: none"> <li>Battery Capacity &gt;10%, NO DC Input</li> <li>No Battery with DC Input</li> </ul> |          |
| Orange Fast Blinking            |    | Power OFF        | ≤10%  | YES      |
| Orange Slow Blinking            |   | Power ON         | ≤10%  | YES      |
| Green Fast Blinking             |  | Power OFF        | >10%  | YES      |
| Green Slow Blinking             |  | Power ON         | >10%  | YES      |
| OFF after Green ON with 1 sec.  |  | Power OFF        | >10%  | NO       |
| OFF after Orange ON with 1 sec. |  | Power OFF        | ≤10%  | NO       |





Link indicator is as table:

| Link Indicator | Lighting Status   | Description  |
|----------------|---|--|
| OFF            |  | <ul style="list-style-type: none"> <li>Panel shut down</li> <li>wired connection broken and wireless connection not ready</li> </ul>   |
| Blue blinking  |  | <ul style="list-style-type: none"> <li>Client mode, wireless connection is ready for connection, but not connected</li> </ul>          |
| Blue ON        |  | <ul style="list-style-type: none"> <li>Client mode, wireless connection is built</li> <li>AP mode, wireless AP is ready for</li> </ul> |




|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 40/137       |

|                |   |  |
|----------------|---|--|
|                |   | connecting   |
| Green ON       |  | <ul style="list-style-type: none"> <li>Wired Connection is built</li> </ul>                            |
| Green blinking |  | <ul style="list-style-type: none"> <li>Panel Initialization</li> <li>Infrared configuration</li> </ul> |

Status indicator is as table:

| Status Indicator | Lighting Status  | Description  |
|------------------|--|--|
| OFF              |   | <ul style="list-style-type: none"> <li>Panel shut down</li> <li>Panel is idle</li> </ul> |
| Green ON         |   | <ul style="list-style-type: none"> <li>Data Transmission</li> </ul>                      |
| Orange blinking  |   | <ul style="list-style-type: none"> <li>Fatal Error</li> </ul>                            |
| Orange ON        |  | <ul style="list-style-type: none"> <li>Panel Initialization</li> </ul>                   |




AP/Client indicator is as table:

| AP/Client Indicator | Lighting Status   | Description  |
|---------------------|---|--|
| OFF                 |  | <ol style="list-style-type: none"> <li>Network connection error;</li> <li>Network connection is wired connection;</li> </ol> |
| Green ON            |  | Wireless connection ok, connection mode is AP  |
| Blue ON             |  | Wireless connection ok, connection mode is Client  |

## 3.2 Battery Charging Installation

|           |        |
|-----------|--------|
| Operating | Figure |
|-----------|--------|

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 41/137       |

|   |  |
|---|--|
| Unload Battery from battery charger.  |    |
| Insert battery into battery charger. Note align the interface position as figure. |   |
| Press the battery to the bottom of battery compartment.                           |  |

### 3.3 Software Installation

In the case of IDemo doesn't work, please install following VC redistribute package.




### 3.4 Panel Infrastructure

Mars1717V supports three connection modes as follows:


- 1) Wired connection mode



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 42/137       |



## 2) Wireless Client Mode



## 3) Wireless AP Mode



To build connection between workstation and Panel, User should follow steps below.

### 3.4.1.1 Wired Mode

To complete Wired connection configuration, users have to finish actions listed below.

#### Configuration of Ethernet Card


To begin configuration of Ethernet Card, users should finish 31.1.2

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 43/137       |


Connect detector to PC with Ethernet Cable





Open Network configuration



Open Local network configuration



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 44/137       |

|                                    |  |
|------------------------------------|--|
| open IPV4 setting                  |    |
| IP setting<br>Network mask setting | IP address: 192.168.8.188<br>Subnet mask: 255.255.255.0                              |
| Open SDK and start connection      |  |
| IP and port setting                | IP: 192.168.8.188<br>Port: 28000   |

### 3.4.1.2 Wireless Client mode


To complete Wireless Client mode configuration, users have to finish actions listed below.

#### Configuration of External wireless AP


|  |   |
|--|---|
| Connect one end of Gigabit Ethernet Cable to PC, Connect another end to LAN port of External wireless AP | / |
|--|---|

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 45/137       |



## Open Network configuration



## Open Local network configuration





|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 46/137       |

|  |  |
|--|--|
| open IPV4 setting  |    |
| IP setting<br>Network mask setting                                 | Select "Obtain an IP address automatically"  |
| Open browser and type 192.168.1.1<br>Log into external wireless AP |  |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 47/137       |

|                                   |  |
|-----------------------------------|--|
| Wireless setup                    | <div> <div>Wireless Setup</div> <div> <div>Apply ▶</div> <div>✕Cancel</div> </div> <hr/> <div>Region Selection</div> <div>Region: Asia ▼</div> <hr/> <div>Wireless Network (2.4GHz b/g/n)</div> <div> <input checked="" type="checkbox"/> Enable SSID Broadcast </div> <div> Name (SSID): NETGEAR_BIG_24 </div> <div> Channel: 13 ▼ </div> <div> Mode: Up to 600 Mbps ▼ </div> <hr/> <div>Security Options</div> <div> <input type="radio"/> None<br/> <input type="radio"/> WPA-PSK [TKIP]<br/> <input checked="" type="radio"/> WPA2-PSK [AES]<br/> <input type="radio"/> WPA-PSK [TKIP] + WPA2-PSK [AES]<br/> <input type="radio"/> WPA/WPA2 Enterprise </div> <hr/> <div>Passphrase: 12345678 (8-63 characters or 64 hex digits)</div> <hr/> <div>Wireless Network (5GHz a/n/ac)</div> <div> <input checked="" type="checkbox"/> Enable SSID Broadcast </div> <div> Name (SSID): NETGEAR_BIG_50 </div> <div> Channel: 153 ▼ </div> <div> Mode: Up to 1300 Mbps ▼ </div> <hr/> <div>Security Options</div> <div> <input type="radio"/> None<br/> <input checked="" type="radio"/> WPA2-PSK [AES]<br/> <input type="radio"/> WPA-PSK [TKIP] + WPA2-PSK [AES]<br/> <input type="radio"/> WPA/WPA2 Enterprise </div> <hr/> <div>Passphrase: 12345678 (8-63 characters or 64 hex digits)</div> <hr/> </div> |
| Configure 2.4GHz wireless network | SSID: NETGEAR_BIG_24<br><br>Security: WPA2-PSK<br><br>Password: 12345678<br><br>Channel: [Please check the current Wi-Fi environment, and choose a relatively clean channel]   |
| Configure 5GHz wireless network   | SSID: NETGEAR_BIG_50<br><br>Security: WPA2-PSK<br><br>Password: 12345678<br><br>Channel: [Please check the current Wi-Fi environment, and choose a relatively clean channel]   |

|  |  |         |            |
|--|--|---------|------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02 |
|  |  | Version | A2         |
|  | Mars1717V User' s Manual               | Date    | 2016.02.23 |
|  |  | Page    | 48/137     |




|   |  |
|---|--|
| LAN setup   |  <p>The screenshot shows the 'LAN Setup' interface. At the top, there are 'Apply' and 'Cancel' buttons. Below is a 'Device Name' field with 'R7000'. The 'LAN TCP/IP Setup' section includes fields for IP Address (192.168.8.1), IP Subnet Mask (255.255.255.0), RIP Direction (Both), and RIP Version (Disabled). There is a checkbox for 'Use Router as DHCP Server' which is checked. Below this are 'Starting IP Address' (192.168.8.2) and 'Ending IP Address' (192.168.8.254). At the bottom is an 'Address Reservation' table with columns for #, IP Address, Device Name, and MAC Address, and buttons for '+ Add', 'Edit', and 'Delete'.</p>   |
| Configure LAN IP address  | IP address: 192.168.8.1<br>Subnet Mask: 255.255.255.0  |
| External Wireless AP Reboot                                       | Apply above settings and reboot your wireless router.  |
| Recover Local Network IPv4 setting of PC wired Ethernet interface |  <p>The screenshot shows the 'Internet Protocol Version 4 (TCP/IPv4) Properties' dialog box. The 'General' tab is selected. It contains instructions about automatic IP assignment. The 'Obtain an IP address automatically' radio button is unselected, and 'Use the following IP address:' is selected. The IP address is 192.168.8.188, Subnet mask is 255.255.255.0, and Default gateway is empty. Below, 'Obtain DNS server address automatically' is unselected, and 'Use the following DNS server addresses:' is selected. Preferred DNS server and Alternate DNS server are empty. There is a checkbox for 'Validate settings upon exit' and an 'Advanced...' button. At the bottom are 'OK' and 'Cancel' buttons.</p> |
| IP setting  | IP address: 192.168.8.188  |
| Network mask setting  | Subnet mask: 255.255.255.0   |

### Configuration of detector

Either Wired Cable or Infrared device can be used to configure detector in wireless client mode.




|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 49/137       |

a.To start configuration with wired cable. It is necessary to finish 3.2.2.1, then proceed to the steps below.

|                            |   |
|----------------------------|---|
| Click “Configure” in IDemo |    |
| Click “Network setting”    |   |
| Click “Add”                |  |



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 50/137       |



|  |  |
|--|--|
| Change SSID and Key,<br>Click “Apply”                          |  |
| Click “write Config to FPD”                                    |  |
| Do not remove wired cable until “write to FPD” button recovers |  |

Since we have chosen default SSID and password, it would connect to wireless AP immediately after powered on next time.



B. To start Infrared configuration. Please see below


|   |   |
|---|---|
| Connect Infrared device with<br>Workstation | / |
|---|---|

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 51/137       |


|  |   |  |
|--|---|--|
| Start IrDARegister.exe   |  |  |
| Click “  ” to open wifi setting | /   |  |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 52/137       |


|  |   |  |
|--|---|--|
| Change SSID and password,do not select AP mode   |  |  |
| Click “Apply”  | /   |  |
| Click”  ” | /   |  |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 53/137       |

|  |   |
|--|---|
| Click"Start"   |  |
| Point Infrared device to detector's infrared interface | /   |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 54/137       |

|                                    |   |
|------------------------------------|---|
| Do not click "Exit" until success  |  |
| Disconnect Infrared device from PC | /   |


### 3.4.1.3 Wireless AP mode

To complete Wired connection configuration, user has to finish actions listed below.






#### Configuration of detector

Either Wired cable or Infrared device can be used to configure detector wireless AP mode.

a. To start wired cable configuration, users should finish 3.4.1.1, then proceed to the steps below.

|                            |   |
|----------------------------|---|
| Click "Configure" in IDemo |  |
|----------------------------|---|


|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 55/137       |

|  |  |
|--|--|
| Click “Network setting”  |   |
| Select AP mode   |   |
| Change SSID and password setting                               |    |
| Click “write Config to FPD”                                    |    |
| Do not remove wired cable until “write to FPD” button recovers | <p>Set FPD WiFi Configuration Done.</p>  |




|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 56/137       |

Since we have chosen default SSID and password, it would connect to wireless AP immediately after powered on next time.

b.To start Infrared configuration, please see below



|                                 |   |
|---------------------------------|---|
| Connect Infrared device with PC | /   |
| Start IrDARegister.exe          |  |

|  |  |         |            |
|--|--|---------|------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02 |
|  |  | Version | A2         |
|  | Mars1717V User' s Manual               | Date    | 2016.02.23 |
|  |  | Page    | 57/137     |



|  |  |
|--|--|
| Click “  ” to open wifi setting |   |
| Select “AP mode”   |  |



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 58/137       |




|  |   |  |
|--|---|--|
| Change SSID and password and other parameter   |  |  |
| Click “Apply”  | /   |  |
| Click”  ” | /   |  |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 59/137       |



|                                  |  |  |
|----------------------------------|--|--|
| Click"Start"                     |   |  |
| Do not click"Exit" until success |  |  |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 60/137       |

### Configuration of external wireless card

|  |  |  |
|--|--|--|
| Open Wireless signal List  |    |  |
| Select SSID which belongs to Detectors, input password and log into system |   |  |
| Open wireless card configuration   |  |  |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 61/137       |

|                               |   |
|-------------------------------|---|
| open IPV4 setting             |   |
| IP setting                    | IP address: 192.168.8.188   |
| Network mask setting          | Subnet mask: 255.255.255.0  |
| Open SDK and start connection |  |
| IP and port setting           | IP: 192.168.8.188   |
|                               | Port: 28000   |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 62/137       |

## 4 Operation

Mars1717V provides SDK for users to integrate detector into their DR system. Additionally, it also provides an application for demonstration, i.e. Idemo. User can use Idemo to control detector without DR system.


### 4.1 Main Operation

To Acquire X-ray image is the main operation of Mars1717V. Most importantly, detector should build synchronization with X-ray generator. Mars1717V is born with two ways to acquire X-ray image, that is Software Mode, Inner Mode and Isync Plus Mode.

#### 4.1.1 Software Mode

##### 4.1.1.1 Block Diagram


Software mode is the basic way to acquire X-ray image. Please see figure below for general feature



Workstation is a host PC device installed with idemo and SDK. Chapter 3.1 has described how to establish connection between detectors and workstation. In software mode, workstation does not have to control X-ray generator. Users would decide when to shoot X-ray.


|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 63/137       |

#### 4.1.1.2 Work flow



#### 4.1.1.3 Timing Setting

To set a clear scenario for program, see diagram below for details



1. Workstation receives "prep" request, send command "Clear" to panel.

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 64/137       |

2. Panel receives “clear” from workstation, start detector internal clear cycle. At the same time, detector would tell workstation “Exposure Prohibited”.
3. Detector finished ”Clear” action and send a message reminding “Exposure Enable”
4. Workstation shows “Exposure Enable” on the iDemo’s message bar to tell user shoot X-ray now.
5. User triggers X-ray generator to initialize and do anode rotation to prepare for X-ray shooting.
6. X-ray generator finishes preparation for X-ray shooting and reminds user to shoot.
7. X-ray generator starts releasing X-ray
8. X-ray generator finishes X-ray shooting.
9. Workstation receives “Acquire” request, send command “Data Acquisition” to panel.
10. Panel receives “Data Acquisition” from workstation, start data acquisition operation.
11. Panel completes image acquisition and begins to send data to workstation.
12. Workstation receives all image data from panel.

If Hardware Pre-offset and Hardware calibration is selected, image is the final image.

If Software Pre-offset and Software Calibration is selected, image would be the raw image, workstation would finish image processing work and show on screen.


If Hardware Post offset and Hardware calibration is selected, image got would be preview image (2x2 binning). After step12, panel would do another dark image acquisition. With both light image and dark image, panel completes all the correction and calibration process. Finally, panel uploads processed image to workstation.

If Software Post offset and Software calibration is selected, image got would be preview image (No binning). After step12, Workstation sends another “clear Acquire” command to panel, panel finishes a dark image acquisition and uploads dark raw image to workstation. With both light image and dark image, workstation completes all the correction and calibration process. Finally, corrected image shows on screen.

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 65/137       |


## 4.1.2 Inner Mode

### 4.1.2.1 Block Diagram



Workstation is a host PC device installed with iDemo and SDK. Chapter 3 has described how to establish connection between panels and workstation. In inner mode, workstation does not control X-ray generator. Users would decide when to shoot X-ray.

### 4.1.2.2 Work Flow






|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 66/137       |

#### 4.1.2.3 Timing Setting

To set a clear scenario for program, see diagram below for details



1. Workstation receives “prep” request and sends “Clear” to panels.
2. Panel receives “clear” from Workstation, start clear operation. Meanwhile, panel would send “Exposure Prohibited” to Workstation.
3. Panel finishes “Clear” operation and send “Exposure Enable” to Workstation.
4. Workstation shows “Exposure Enable” on the iDemo’s message bar to tell user shoot X-ray.
5. User triggers X-ray generator to initialize and do anode rotation to prepare for X-ray shooting
6. X-ray generator finishes preparation and reminds users.
7. X-ray generator begins releasing X-ray
8. X-ray generator finishes X-ray shooting.
9. X-ray sensor in panel triggers panel to start image acquisition operation.
10. Panel completes image acquisition and begins to send data to Workstation.
11. Workstation starts receiving image data from panel.
12. Workstation receives all image data from panel.

If Hardware Pre-offset and Hardware calibration is selected, image got is the final image.

If Software Pre-offet and Software Calibration is selected, image got would be raw image, workstation would finish image processing work and show on screen.

If Hardware Post offset and Hardware calibration is selected, image got from detector would be preview image (2x2 binning). After step12, Detector would do another dark image acquisition. With both light image

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 67/137       |

and dark image, detector completes all the correction process. Finally, detector uploads corrected image and workstation shows on screen.

If Software Post offset and Software calibration is selected, image got from panel would be preview image (No binning). After step12, Workstation sends another “clear Acquire” to panel, panel would do dark image acquisition and uploads dark raw image to workstation. With both light image and dark image, workstation completes all the correction process. Finally, corrected image shows on screen.


#### 4.1.2.4 Abnormal Action

Action1: after Step4, if user wants to cancel this exposure cycle, Idemo provides an “Abort Exp” function to close exposure windows. However, Idemo allows user to click “Abort Exp” until Workstation receives first image.

Action2: after Step4, if user does not shoot X-ray in exposure windows, panel would close exposure windows automatically and send a message to workstation that waiting for X-ray shooting is overtime. Meanwhile, panel would also start image acquisition. After image acquisition, panel sends image to workstation.

### 4.1.3 Isync Plus Mode

#### 4.1.3.1 Block Diagram




|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 68/137       |


Workstation is a host PC device installed with iDemo and SDK. Chapter 3 has described how to establish connection between panel and Workstation. In Isync Plus mode, User doesn't interact with Workstation.

After shooting, images would be shown on screen immediately.

#### 4.1.3.2 Work Flow



#### 4.1.3.3 Timing Setting



1. X-ray generator is ready for X-ray shooting and begins to release X-ray.

2. Workstation receives "Exposure Prohibited" from Panel.



3. Panel starts uploading Pre-dark image and Light image to Workstation for preview. If hardware offset is selected, panel would do offset first, and then upload preview image (2X2 binning).

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 69/137       |

4. Panel starts uploading Post-dark image to Workstation. If hardware offset is chosen, panel would do correction and calibration first, then upload processed image to Workstation.

5. Workstation receives “Exposure Enable” from Panel.

## 4.2 Connection Build

|   |  |
|---|--|
| Click “Start”   |    |
| Input IP address and port number. The IP address should be the same as the IP address of the network card connected with panel. The port should use the default value of 28000;<br>Click “Connect”; |  |



Note:

1. Once changing connection from different network card, user must re-connect panel with different IP address.
2. Switching between wired and wireless connection does not need any extra operation.

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 70/137       |

3. The rule of Multi-Share control is based on IP address. The second terminal with different IP address is not allowed to operate panel after the first one connected. If there is no command transmission between panel and Workstation over 5 minutes, panel releases access authority.

## 4.3 Panel Configuration

|                |   |
|----------------|---|
| Click "Config" |    |
| Set parameters |  |

Note: 1. If panel works in Isync Plus mode, it is not allowed to change any parameters and write into ROM or RAM. User is required to switch to software mode, change parameters and then switch to Isync Plus mode. On the other hand, we do not recommend user to switch working mode too often.


## 4.4 Correction and Calibration Template Generation

Iray recommends performing correction and calibration after installation or any major change on the system settings and hardware configuration. On the other hand, it is also recommended to do the correction and calibration in each 6 months.

### 4.4.1 Pre-offset Template Generation


If panel is configured to do Pre-offset correction, Pre-offset Template is necessary. See below

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 71/137       |




|  |  |
|--|--|
| Click “ Create”, Choose “Defect”                         |  |
| Click “Create Offset”, wait until image acquisition ends | /  |

#### 4.4.2 Gain Calibration Template Generation


Before Gain template generating, make sure SID1.2m, no copper is required,

|  |  |
|--|--|
| Choose “Post” offset mode <sup>1</sup> |  |
|--|--|

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 72/137       |

|   |  |
|---|--|
| Click “ Create”, Choose “Gain”  |    |
| Set “Expected Gray Value” 10000;<br>Click “New Gain Point”;<br>Set “Acquire Delay” 3000;<br>If panel is in software mode, click “Prep”, shoot X-ray, click “Acquire”;<br>If panel is in Inner mode, click “Prep”, shoot X-ray;<br>If panel is in Isync Plus mode, click “Prep Acquire”, wait for ready, shoot X-ray. Time window can be set by “Acquire Delay”; |   |
| Wait for Post-offset image uploaded, there will be value shown in “current average” box, change X-ray dose to make sure “current average” in the range of $10000 \pm 100$ , If yes, Click “Select”  |  |


|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 73/137       |

|  |  |
|--|--|
| <p>Repeat the same operation for 4 times, “Create Gain Correct file” button can be available, Click “Create Gain Correct file”</p> |  |
| <p>Click “OK”</p>  | <p>/</p>   |

Note:1 please use software post offset correction.

### 4.4.3 Defect Correction Template Generation


Before Defect template generating, make sure SID1.2m, no copper is required,

|  |  |
|--|--|
| <p>Choose “Post” offset mode<sup>1</sup></p> |  |
|--|--|




|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 74/137       |


Click “ Create”, Choose “Defect”



Click “start creating”, message box will show you the first image value for defect correction<sup>2</sup>;  
If panel is in software mode, click “Prep”, shoot X-ray, click “Acquire”;  
If panel is in Inner mode, click “Prep”, shoot X-ray;  
If panel is in Isync Plus mode, click “Prep Acquire”, wait for ready, shoot X-ray. Time window can be set by “Acquire Delay”;




Click “start creating” to start the second X-ray shoot




|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 75/137       |


Click “start creating” to start the third X-ray shoot




Click “start creating” to start the fourth X-ray shoot



After the fourth X-ray shoot, you do not need to change the X-ray dose, just repeat operation of “start creating” and image acquisition until it comes to the 19th images.



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 76/137       |


|   |  |
|---|--|
| Click “Create Defect”, wait until it ends |  |
|---|--|

Note:1 please use software post offset mode.



2. Make sure your X-ray dose is right, if your dose is out of the range, idemo will remind you to adjust the dose. Then you can click “start creating” and try again.

#### 4.4.4 Lag Template Generation

Lag template is only used in Isync plus mode, if panel works in other mode, user does not have to complete lag template. Before Lag template generating, make sure SID 1.2m, no copper is required.

|   |  |
|---|--|
| Choose “Post” offset mode<br>Choose “curr” gain mode<br>Choose “curr” defect mode | /  |
| Click “Create”<br>Choose “Lag” correction<br>Click “Start creating”               |  |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 77/137       |

|   |   |
|---|---|
| <p>Change X-ray dose to make sure average gray value in the range of <math>10000 \pm 100</math>, Shoot X-ray, Post-offset image is shown on screen, make sure exposure line between 100 and 500;</p> <p>Shoot X-ray, another Post-offset image is shown on screen, make sure exposure line between 2572 and 2972;</p> |   |
| <p>Image acquisition finished; Click "Create Lag Map";</p>  |  |
| <p>Display "Lag calibration finished", correction finished, close correction window</p>   | /   |

## 4.5 Image Acquisition Continually


Image Acquisition continually can be used only in software mode and inner mode, it is not supported in Isync Plus mode. The operation is designed for panel testing, not for customer using.

"SEQU" is the command to start image acquisition and "STOP" is the opposite. If user wants to change frequency of image acquisition, change parameters in "Config-General Setting-Time Span".


## 4.6 Image save

"SAVE" provides two features of image saving for user. The first is saving multiple images, the second is saving single image.



### 4.6.1 Multiple images

|              |   |
|--------------|---|
| Click "Save" |  |
|--------------|---|


|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 78/137       |

|   |  |
|---|--|
| <p>Set the number of images in blank;<br/>Choose the saving path;<br/>Click “Begin Save”;</p> |  |
|---|--|

#### 4.6.2 Single image

|                              |  |
|------------------------------|--|
| <p>Click “Save”</p>          |    |
| <p>Click “Save Directly”</p> |  |



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 79/137       |

|  |  |
|--|--|
| <p>Choose the saving path;<br/>Click “save”;</p> |  |
|--|--|

## 4.7 Image Check and upload

“OPEN” provides three feature for image check and uploading. Local Image Check, Panel Image Upload and Panel Image Check. Local Image Check defines function to check image saved in Workstation. Panel Image Upload defines function to upload images stored in panel. Panel Image Check defines function to check images stored in panel.

### 4.7.1 Local Image Check


|   |  |
|---|--|
| <p>Click “Open”</p>   |   |
| <p>Click “Open Local File”</p>  |  |
| <p>Choose images stored in Workstation, images would be shown on screen</p> | /  |

|  |  |                |              |
|--|--|----------------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | <b>Doc NO.</b> | 037-201-02   |
|  |  | <b>Version</b> | A2           |
|  | Mars1717V User' s Manual               | <b>Date</b>    | 2016. 02. 23 |
|  |  | <b>Page</b>    | 80/137       |

## 4.7.2 Panel Image Upload

Before panel image uploading, FTP server is necessary.

Make sure firewall is closed




Start "FTPServer.exe"

| 名称                  | 修改日期             | 类型     | 大小       |
|---------------------|------------------|--------|----------|
| rpaCamLinkAcqu.dll  | 2015/7/6 18:20   | 应用程序扩展 | 8 KB     |
| FpdCommonFunc.dll   | 2015/7/6 18:20   | 应用程序扩展 | 31 KB    |
| FpdCorrection.dll   | 2015/7/6 18:20   | 应用程序扩展 | 53 KB    |
| FpdDefect.dll       | 2015/7/6 18:20   | 应用程序扩展 | 44 KB    |
| FpdDetector.dll     | 2015/7/6 18:20   | 应用程序扩展 | 10 KB    |
| FpdDetectorE3.dll   | 2015/7/6 18:20   | 应用程序扩展 | 38 KB    |
| FpdDetectorE4.dll   | 2015/7/6 18:20   | 应用程序扩展 | 71 KB    |
| FpdFTS.dll          | 2015/7/6 18:20   | 应用程序扩展 | 8 KB     |
| FpdImageProcess.dll | 2015/7/6 18:21   | 应用程序扩展 | 57 KB    |
| FpdIrDARegister.dll | 2015/5/21 8:59   | 应用程序扩展 | 16 KB    |
| FpdKey.dll          | 2015/7/6 18:20   | 应用程序扩展 | 8 KB     |
| FpdLib.dll          | 2015/7/6 18:21   | 应用程序扩展 | 83 KB    |
| FpdSyncbox.dll      | 2015/7/6 18:20   | 应用程序扩展 | 16 KB    |
| FpdTCP.dll          | 2015/7/6 18:20   | 应用程序扩展 | 12 KB    |
| FpdUDP.dll          | 2015/7/6 18:20   | 应用程序扩展 | 13 KB    |
| FTPServer.exe       | 2014/5/20 0:10   | 应用程序   | 123 KB   |
| GridLineRemoval.dll | 2014/2/7 10:38   | 应用程序扩展 | 219 KB   |
| IDemo_interface.exe | 2015/7/6 18:21   | 应用程序   | 1,006 KB |
| IDicom.dll          | 2013/11/18 13:44 | 应用程序扩展 | 2,976 KB |
| IrDARegister.exe    | 2015/6/9 9:42    | 应用程序   | 130 KB   |
| log4cpp32vc9.dll    | 2014/11/11 10:19 | 应用程序扩展 | 237 KB   |
| users.dat           | 2014/5/30 9:17   | DAT 文件 | 1 KB     |


|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 81/137       |

Choose “Configuration”;  
Set FTP Port “21000”;  
Click “Start” button;




Panel Image is uploaded as following.

Click “Open”




Click “Get File List”






|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 82/137       |


Images stored in panel would be shown in the list. From the list, user could check basic information of images



Select images user wants;  
Click “Start Upload”;  
Images would be uploaded to  
\*/idemo/upload





If user want to upload multiple images,  
select multiple images, click “start  
upload”, during uploading, user could  
click “stop upload” to stop uploading



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 83/137       |

### 4.7.3 Panel Image Check


If user wants to check images stored in panel immediately, see below

|   |   |
|---|---|
| Click "Open"                                    |    |
| Select image user need;<br>Click "open in list" |  |




## 4.8 Defect Template Check and Modification

Idemo provides function to check defect template. If defect template has updates, user could add and delete defect pixel or defect lines by modifying defect template opened.


### 4.8.1 Defect Template Check

|              |   |
|--------------|---|
| Click "Open" |  |
|--------------|---|

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 84/137       |


|   |   |
|---|---|
| Select defect template;<br>Click “open”;  |   |
| Select Product Size “Venu1417 (2304*2800);<br>Click “OK”, Defect template is shown on screen; |  |
| Click “Close” to close Defect template  |  |

## 4.8.2 Defect Template Modification


|              |   |
|--------------|---|
| Click “Open” |  |
|--------------|---|

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 85/137       |


Select Defect Template;  
Click “Open”;





Set product size “Venu1417 (2304\*2800);  
Click “OK”



If there is new defect pixel, input  
coordinate, click “Add”;  
If pixel is labeled as defect by mistake,  
input coordinate, click “Delete”




|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 86/137       |

|   |  |
|---|--|
| <p>If there is new defect line, input coordinate, click “Add”;</p> <p>If line is labeled as defect by mistake, input coordinate, click “Delete”</p> |  |
| Click “Close”   |   |



## 4.9 Correction and Calibration Management


### 4.9.1 Correction and Calibration template synchronization

Panel supports correction and calibration template storage. So template in panel could be uploaded to Workstation, and template in Workstation could also be downloaded to panel. Before synchronization, FTP server is necessary.



|                              |  |
|------------------------------|--|
| Make sure firewall is closed |  |
|------------------------------|--|

|  |  |                |              |
|--|--|----------------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | <b>Doc NO.</b> | 037-201-02   |
|  |  | <b>Version</b> | A2           |
|  | Mars1717V User' s Manual               | <b>Date</b>    | 2016. 02. 23 |
|  |  | <b>Page</b>    | 87/137       |

|   |   |
|---|---|
| Start “FTPServer.exe”   |   |
| Choose “Configuration”;<br>Set FTP Port “21000”;<br>Click “Start” button; |  |

|                |   |
|----------------|---|
| Click “Config” |  |
|----------------|---|


|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 88/137       |

|  |   |
|--|---|
| Select “HW Correction”   |   |
| <p>If user wants to download correction and calibration template to panel. Click “Download” from top to bottom;</p> <p>If user wants to upload correction and calibration template from panel. Click “Read out” and “Update Local” from top to bottom;</p> |  |




## 4.9.2 Correction and Calibration management

Panel supports two ways to do correction and calibration. Software Correction and Calibration defines the scenario that Workstation completes all correction and calibration. If panel complete all correction and calibration by itself, it is named as Hardware Correction and Calibration.


### 4.9.2.1 Software Correction and Calibration

|                |   |
|----------------|---|
| Click “Config” |  |
|----------------|---|

|  |  |                |              |
|--|--|----------------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | <b>Doc NO.</b> | 037-201-02   |
|  |  | <b>Version</b> | A2           |
|  | Mars1717V User' s Manual               | <b>Date</b>    | 2016. 02. 23 |
|  |  | <b>Page</b>    | 89/137       |




|  |  |
|--|--|
| Select “HW Correction”;                                      |    |
| Set Offset Mode “NO”;  |   |
| Set Offset “POST”;<br>Set Gain “CURR”;<br>Set Defect “CURR”; |  |

#### 4.9.2.2 Hardware Correction and Calibration Management

|                |   |
|----------------|---|
| Click “Config” |  |
|----------------|---|



|  |  |                |              |
|--|--|----------------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | <b>Doc NO.</b> | 037-201-02   |
|  |  | <b>Version</b> | A2           |
|  | Mars1717V User' s Manual               | <b>Date</b>    | 2016. 02. 23 |
|  |  | <b>Page</b>    | 90/137       |

|  |  |
|--|--|
| <p>Select “HW Correction”;</p>                 |    |
| <p>Set Offset Mode “Post offset”;</p>          |   |
| <p>Set Gain “CURR”;<br/>Set Defect “CURR”;</p> |  |



## 4.10 Sleep and Wake Up

Panel supports sleep and wake up operation. User can trigger sleep manually or automatically, but there is only one way to wake up panel.

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 91/137       |


## 4.10.1 Sleep

### 4.10.1.1 Manual Sleep




|   |   |
|---|---|
| Click “Sleep”, panel will go to sleep immediately |    |
| Message box shows that panel is sleeping          |  |

### 4.10.1.2 Automatic Sleep

To go to sleep automatically, user should set the time flag first. If panel detects that there is no operation in time flag, panel would go to sleep. If time flag is set zero, panel would not go to sleep automatically.



|                |   |
|----------------|---|
| Click “Config” |  |
|----------------|---|

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 92/137       |

|  |  |
|--|--|
| Set “auto sleep timeout” value user need                                 |    |
| If there is no operation in time flag user sets, panel would go to sleep |   |
| Message box shows that panel is sleeping                                 |  |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 93/137       |

### 4.10.2 Wake Up

|                                       |   |
|---------------------------------------|---|
| If panel is sleeping, click “Wake Up” |    |
| Message box shows that panel is ready |  |

## 4.11 Firmware Update

Panel supports updating firmware. If user wants to update new firmware, see below


### 4.11.1 MCU Update

If current MCU version is 2.5.1.\*, we should follow instruction below.

|  |   |
|--|---|
| Open “mini Cover”  | / |
| Remove original Ethernet cable, insert J-link download cable | / |
| Start J-flash ARM  | / |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 94/137       |


Click “file->open data file”



Click “option->project setting”,  
Set “connection to J-link” USB mode




Click “Target Interface”, Choose “SWD”



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 95/137       |

Click “CPU”;  
Select “ST STM32L151C8”;



Click “Target->Connect”

/

Click “Target-> Erase”

/

Click “Target->Program”

/

Click “Target->Start Application”


/

Note:1. Make sure panel is powered up.




If current MCU version is 2.5.2.\*, Please refer to 4.10.2 and 4.10.3 for upgrading.

#### 4.11.2 FTP Server

Make sure firewall is closed





|  |  |                |            |
|--|--|----------------|------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | <b>Doc NO.</b> | 037-201-02 |
|  |  | <b>Version</b> | A2         |
|  | Mars1717V User' s Manual               | <b>Date</b>    | 2016.02.23 |
|  |  | <b>Page</b>    | 96/137     |

| Click On “FTPServer.exe”  |  <table><thead><tr><th>名称</th><th>修改日期</th><th>类型</th><th>大小</th></tr></thead><tbody><tr><td>FpdCommonFunc.dll</td><td>2015/6/16 10:49</td><td>应用程序扩展</td><td>30 KB</td></tr><tr><td>FpdCorrection.dll</td><td>2015/6/16 10:50</td><td>应用程序扩展</td><td>53 KB</td></tr><tr><td>FpdDefect.dll</td><td>2015/6/16 10:49</td><td>应用程序扩展</td><td>44 KB</td></tr><tr><td>FpdDetector.dll</td><td>2015/6/16 10:49</td><td>应用程序扩展</td><td>10 KB</td></tr><tr><td>FpdDetectorE3.dll</td><td>2015/6/16 10:50</td><td>应用程序扩展</td><td>39 KB</td></tr><tr><td>FpdDetectorE4.dll</td><td>2015/6/16 10:50</td><td>应用程序扩展</td><td>71 KB</td></tr><tr><td>FpdFTS.dll</td><td>2015/6/16 10:50</td><td>应用程序扩展</td><td>8 KB</td></tr><tr><td>FpdImageProcess.dll</td><td>2015/6/16 10:50</td><td>应用程序扩展</td><td>58 KB</td></tr><tr><td>FpdIrDARegister.dll</td><td>2015/5/21 8:59</td><td>应用程序扩展</td><td>16 KB</td></tr><tr><td>FpdKey.dll</td><td>2015/6/16 10:50</td><td>应用程序扩展</td><td>8 KB</td></tr><tr><td>FpdLib.dll</td><td>2015/6/16 10:50</td><td>应用程序扩展</td><td>82 KB</td></tr><tr><td>FpdSyncbox.dll</td><td>2015/6/16 10:50</td><td>应用程序扩展</td><td>16 KB</td></tr><tr><td>FpdTCP.dll</td><td>2015/6/16 10:49</td><td>应用程序扩展</td><td>13 KB</td></tr><tr><td>FpdUDP.dll</td><td>2015/6/16 10:49</td><td>应用程序扩展</td><td>13 KB</td></tr><tr><td>FTPServer.exe</td><td>2014/5/20 0:10</td><td>应用程序</td><td>123 KB</td></tr><tr><td>GridLineRemoval.dll</td><td>2014/2/7 10:38</td><td>应用程序扩展</td><td>219 KB</td></tr><tr><td>IDemo_interface.exe</td><td>2015/6/16 10:52</td><td>应用程序</td><td>1,000 KB</td></tr><tr><td>IDicom.dll</td><td>2013/11/18 13:44</td><td>应用程序扩展</td><td>2,976 KB</td></tr><tr><td>IrDARegister.exe</td><td>2015/6/9 9:42</td><td>应用程序</td><td>130 KB</td></tr><tr><td>log4cpp32vc9.dll</td><td>2014/11/11 10:19</td><td>应用程序扩展</td><td>237 KB</td></tr><tr><td>users.dat</td><td>2014/5/30 9:17</td><td>DAT 文件</td><td>1 KB</td></tr></tbody></table> | 名称     | 修改日期     | 类型 | 大小 | FpdCommonFunc.dll | 2015/6/16 10:49 | 应用程序扩展 | 30 KB | FpdCorrection.dll | 2015/6/16 10:50 | 应用程序扩展 | 53 KB | FpdDefect.dll | 2015/6/16 10:49 | 应用程序扩展 | 44 KB | FpdDetector.dll | 2015/6/16 10:49 | 应用程序扩展 | 10 KB | FpdDetectorE3.dll | 2015/6/16 10:50 | 应用程序扩展 | 39 KB | FpdDetectorE4.dll | 2015/6/16 10:50 | 应用程序扩展 | 71 KB | FpdFTS.dll | 2015/6/16 10:50 | 应用程序扩展 | 8 KB | FpdImageProcess.dll | 2015/6/16 10:50 | 应用程序扩展 | 58 KB | FpdIrDARegister.dll | 2015/5/21 8:59 | 应用程序扩展 | 16 KB | FpdKey.dll | 2015/6/16 10:50 | 应用程序扩展 | 8 KB | FpdLib.dll | 2015/6/16 10:50 | 应用程序扩展 | 82 KB | FpdSyncbox.dll | 2015/6/16 10:50 | 应用程序扩展 | 16 KB | FpdTCP.dll | 2015/6/16 10:49 | 应用程序扩展 | 13 KB | FpdUDP.dll | 2015/6/16 10:49 | 应用程序扩展 | 13 KB | FTPServer.exe | 2014/5/20 0:10 | 应用程序 | 123 KB | GridLineRemoval.dll | 2014/2/7 10:38 | 应用程序扩展 | 219 KB | IDemo_interface.exe | 2015/6/16 10:52 | 应用程序 | 1,000 KB | IDicom.dll | 2013/11/18 13:44 | 应用程序扩展 | 2,976 KB | IrDARegister.exe | 2015/6/9 9:42 | 应用程序 | 130 KB | log4cpp32vc9.dll | 2014/11/11 10:19 | 应用程序扩展 | 237 KB | users.dat | 2014/5/30 9:17 | DAT 文件 | 1 KB |
|---|---|--------|----------|----|----|-------------------|-----------------|--------|-------|-------------------|-----------------|--------|-------|---------------|-----------------|--------|-------|-----------------|-----------------|--------|-------|-------------------|-----------------|--------|-------|-------------------|-----------------|--------|-------|------------|-----------------|--------|------|---------------------|-----------------|--------|-------|---------------------|----------------|--------|-------|------------|-----------------|--------|------|------------|-----------------|--------|-------|----------------|-----------------|--------|-------|------------|-----------------|--------|-------|------------|-----------------|--------|-------|---------------|----------------|------|--------|---------------------|----------------|--------|--------|---------------------|-----------------|------|----------|------------|------------------|--------|----------|------------------|---------------|------|--------|------------------|------------------|--------|--------|-----------|----------------|--------|------|
| 名称  | 修改日期  | 类型     | 大小       |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| FpdCommonFunc.dll   | 2015/6/16 10:49   | 应用程序扩展 | 30 KB    |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| FpdCorrection.dll   | 2015/6/16 10:50   | 应用程序扩展 | 53 KB    |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| FpdDefect.dll   | 2015/6/16 10:49   | 应用程序扩展 | 44 KB    |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| FpdDetector.dll   | 2015/6/16 10:49   | 应用程序扩展 | 10 KB    |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| FpdDetectorE3.dll   | 2015/6/16 10:50   | 应用程序扩展 | 39 KB    |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| FpdDetectorE4.dll   | 2015/6/16 10:50   | 应用程序扩展 | 71 KB    |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| FpdFTS.dll  | 2015/6/16 10:50   | 应用程序扩展 | 8 KB     |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| FpdImageProcess.dll   | 2015/6/16 10:50   | 应用程序扩展 | 58 KB    |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| FpdIrDARegister.dll   | 2015/5/21 8:59  | 应用程序扩展 | 16 KB    |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| FpdKey.dll  | 2015/6/16 10:50   | 应用程序扩展 | 8 KB     |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| FpdLib.dll  | 2015/6/16 10:50   | 应用程序扩展 | 82 KB    |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| FpdSyncbox.dll  | 2015/6/16 10:50   | 应用程序扩展 | 16 KB    |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| FpdTCP.dll  | 2015/6/16 10:49   | 应用程序扩展 | 13 KB    |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| FpdUDP.dll  | 2015/6/16 10:49   | 应用程序扩展 | 13 KB    |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| FTPServer.exe   | 2014/5/20 0:10  | 应用程序   | 123 KB   |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| GridLineRemoval.dll   | 2014/2/7 10:38  | 应用程序扩展 | 219 KB   |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| IDemo_interface.exe   | 2015/6/16 10:52   | 应用程序   | 1,000 KB |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| IDicom.dll  | 2013/11/18 13:44  | 应用程序扩展 | 2,976 KB |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| IrDARegister.exe  | 2015/6/9 9:42   | 应用程序   | 130 KB   |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| log4cpp32vc9.dll  | 2014/11/11 10:19  | 应用程序扩展 | 237 KB   |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| users.dat   | 2014/5/30 9:17  | DAT 文件 | 1 KB     |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| Choose “Configuration”, set “FTP Port” 21000, others as default |  <p>Quick'n Easy FTP Server Lite</p> <p>Server View Help</p> <p>Configuration</p> <p>FTP Port: 21000</p> <p>Max Connections: 10</p> <p>Connection Timeout: 5 minutes</p> <p>Welcome Message: Welcome to Pablo's FTP Server</p> <p>Goodbye Message: Bye</p> <p><input type="checkbox"/> Launch FTP Server at windows startup</p> <p><input type="checkbox"/> Startup minimized in systemtray</p> <p><input type="checkbox"/> Automatically activate server at startup</p> <p>Apply</p> <p>FTP Server is offline 0 bytes received 0 bytes sent</p>   |        |          |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |
| Start FTP Server  |  <p>Quick'n Easy FTP Server Lite</p> <p>Server View Help</p> <p>Server Log</p> <p>FTP Server failed to listen on port 21000.</p> <p>FTP Server is offline 0 bytes received 0 bytes sent</p>   |        |          |    |    |                   |                 |        |       |                   |                 |        |       |               |                 |        |       |                 |                 |        |       |                   |                 |        |       |                   |                 |        |       |            |                 |        |      |                     |                 |        |       |                     |                |        |       |            |                 |        |      |            |                 |        |       |                |                 |        |       |            |                 |        |       |            |                 |        |       |               |                |      |        |                     |                |        |        |                     |                 |      |          |            |                  |        |          |                  |               |      |        |                  |                  |        |        |           |                |        |      |



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 97/137       |

### 4.11.3 Firmware Update

|   |   |
|---|---|
| Put update file in the directory of<br>“*\\idemo\\update” |   |
| Click “Config”, Choose Firmware Update                    |  |



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 98/137       |

|   |   |
|---|---|
| Choose Update Type “core” <sup>1</sup>  |   |
| Click “Select”, choose right update file <sup>2</sup>   |   |
| Click “update”, waiting for message box <sup>3</sup>  |  |
| Click message box, waiting for end of rebooting <sup>4</sup> ;<br>Panel finishes rebooting, Click “OK” <sup>5</sup> | /   |

Note:

1. It is not limited to “Core”, actually, other choice is also ok.
2. If it is MCU update, choose MCU image file. Otherwise, choose ALL-Image file, Please make sure update file is selected, if not, panel will be not in use after updating.
3. There is a progress bar for indication. Make sure battery is inserted and battery capacity is over 25%
4. This rebooting function is controlled by panel itself. It has the same function with “ Reset FPD”
5. Please make sure Idemo show “ Ready”. It can also be checked by click “ Config” button, there is firmware version.

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 99/137       |

## 4.12Short cut


iDemo supports some shortcuts as follows:

- Double-click the left mouse button, the image displayed in center and with maximum size.
- Double-click the right mouse button, the window level and width adjusted to WL:8191/WW:16383.
- Drag the left mouse button, drag the image displayed.
- Lateral-drag the right mouse button to adjust the window width, and vertical-drag the right mouse button to adjust the window level.
- F3 Key: Quickly adjust the image window width and window level.

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 100/137      |

## 4.13 Software

### 4.13.1 Main GUI



iRay provides test tools, such as iDemo for testing the basic performance of detector. It can connect the detector, acquire image, image correct and calibrate.

Function description of regions and buttons within the main window as follows:

|   |            |                                    |    |            |  |
|---|------------|------------------------------------|----|------------|--|
| 1 | Start/End  | Load or unload NIC device driver   | 13 | Open/Close | Open or close defect map                   |
| 2 | Sleep/Wake | Sleep or wake up panel             | 14 | Offset     | Open or close software post offset         |
| 3 | Prep       | Clear lags of the panel            | 15 | Gain       | Open or close gain calibration             |
| 4 | Abort Exp  | Close exposure Window              | 16 | Defect     | Open or close defect correction            |
| 5 | Acquire    | Acquire an image without clearance | 17 | Create     | Generate gain template and defect template |

|  |  |         |            |
|--|--|---------|------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02 |
|  |  | Version | A2         |
|  | Mars1717V User' s Manual               | Date    | 2016.02.23 |
|  |  | Page    | 101/137    |

|    |              |   |    |  |
|----|--------------|---|----|--|
| 6  | Check ID     | Check panel SN                                      | 18 | Version of the idemo                                       |
| 7  | Prep Acquire | Clear lags and acquire an image                     | 19 | Status of the idemo  |
| 8  | Config       | Configure the panel                                 | 20 | Panel feedback message                                     |
| 9  | Sequ         | Start acquiring images continually                  | 21 | Acquisition interval between two images                    |
| 10 | Stop         | Stop acquiring images continually                   | 22 | Pixel X/Y coordinate and gray scale value of Pixel(14 bit) |
| 11 | Save         | Save images continually or save an image            | 23 | Image acquisition instruction box                          |
| 12 | Open         | Open local images(.DCM) or open images in the panel | 24 | WiFi signal and battery capacity indication                |
|    |              |   | 25 | Region of image display                                    |

### 4.13.2 Message Box

#### 4.13.2.1 Status Box

Status box defines the current status of panel.

| Value    | Description   |
|----------|---|
| Offline  | Idemo loose connection with panel, it does not receive heart beat                                     |
| Ready    | Idemo builds connection with panel, panel is ready for receiving new operation                        |
| Busy     | Idemo builds connection with panel, panel is busying on the last operation, it can not be interrupted |
| Sleeping | Idemo builds connection with panel, panel has gone in sleep   |
| Waking   | Idemo builds connection with panel, panel is being wake up  |
| Timeout  | Command executes overtime   |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 102/137      |

#### 4.13.2.2 Feedback Box

Feedback box shows feedback message from panel.

#### 4.13.2.3 Acquisition Interval Box

Acquisition Interval Box shows the time between two image acquired currently.

#### 4.13.2.4 Coordinate and Gray Scale Box

Coordinate and Cray Scale Box show the coordinate and gray scale of mouse.

#### 4.13.2.5 Image Acquisition Box

Image Acquisition Box shows whether image is uploading

| Color | Description                   |
|-------|-------------------------------|
| Red   | Image is uploading from panel |
| Green | others                        |


#### 4.13.2.6 Battery and Connection Box

Battery and Connection Box shows battery capacity, wireless signal level and wire connection.

#### 4.13.2.7 Progress Bar

Progress Bar defines as following.


|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 103/137      |



If progress bar is Green when shooting X-ray, image quality is acceptable, otherwise image quality would degrade.

### 4.13.3 Configuration GUI

#### 4.13.3.1 General Settings



|                | Item           | Description                       | Modify |
|----------------|----------------|-----------------------------------|--------|
| Detector Info  | Product type   | Type of panel product             | NO     |
|                | Serial Number  | Serial number of the panel        | NO     |
| Version Detail | Lib Version    | Version number of idemo           | NO     |
|                | Core Version   | Version number of ARM application | NO     |
|                | Kernel Version | Version number of ARM OS Kernel   | NO     |
|                | Main Version   | Version number of Core FPGA       | NO     |


|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 104/137      |

|                  |                    |   |     |
|------------------|--------------------|---|-----|
|                  | Read Version       | Version number of Read FPGA   | NO  |
|                  | MCU Version        | Version number of MCU   | NO  |
| Sensor details   | Temperature        | Panel inner temperature(Read Board and Core Board)  | NO  |
|                  | Humidity           | Panel inner Humidity  | NO  |
| General Settings | Trigger            | Detector trigger mode:<br>1.Outer<br>2.Inner<br>3.Software<br>4.PREP<br>5.Service<br>6.FreeSync (Default)                       | YES |
|                  | Time Span          | This parameter is used only in continual acquisition. The time span is the time interval between two nearby acquisition process | YES |
|                  | Exp Window         | Exposure window is used in Inner mode, it defines the time for X-ray shooting.  | YES |
|                  | FreeSync Window    | Not Used  | YES |
|                  | Acquire Delay      | Delay time before image acquisition.  | YES |
|                  | Delay time         | Exposure Window in Isync Plus mode/Delay time between clearance and acquisition in other mode.                                  | YES |
|                  | Auto Sleep Timeout | Time span of idle before going to sleep   | YES |
|                  | Integrate          | The integration time for the photo diode  | NO  |
|                  | Inner Delay        | The real delay time between clearance and acquisition most recently.  | NO  |
|                  |                    |   |     |
| Button           | Recover ROM        | Recover configuration to factory setting  | /   |
|                  | Write ROM          | Write configuration data into nonvolatile memory  | /   |
|                  | Read ROM           | Read configuration data from nonvolatile memory   | /   |
|                  | Write RAM          | Write configuration data into volatile memory   | /   |
|                  | Read RAM           | Read configuration data into volatile memory  | /   |
|                  | Reset FPD          | Reboot the panel  | /   |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 105/137      |

|  |                 |                                |   |
|--|-----------------|--------------------------------|---|
|  | Exit            | Exit configuration GUI         | / |
|  | Query Live Time | Check the active time of panel | / |

#### 4.13.3.2 Advanced Settings



|                   | Item                       | Description                                      | Modify |
|-------------------|----------------------------|--|--------|
| Advanced Settings | Row Pre Delay              | Delay time before acquiring row data             | YES    |
|                   | VT                         | Voltage corresponding to the charge compensation | YES    |
|                   | Freesync Center Threshold  | Not used   | YES    |
|                   | Row Post Delay             | Delay time after acquiring row data              | YES    |
|                   | PGA                        | Integrator capacitor range.                      | YES    |
|                   | Freesync Edge Threshold    | Not used   | YES    |
|                   | Detector Auto Clear        | Set the detector in auto clear mode              | NO     |
|                   | Detector Auto Clear Period | Auto clear period for panel                      | NO     |
| Mammo             | Exposure Mode              | Not used   | NO     |
|                   | AEC Main Time              | Not used   | NO     |




|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 106/137      |

|              |                  |  |     |
|--------------|------------------|--|-----|
|              | Dyna Offset Mode | Not used                               | NO  |
|              | Dyna Offset Time | Not used                               | NO  |
|              | Correct Folder   | Not used                               | NO  |
| Shock Sensor | Read Grid Status | Not used                               | NO  |
|              | Read Shock log   | Read shock sensor log                  | NO  |
|              | Clear Shock Log  | Clear shock sensor log                 | NO  |
|              | Set Threshold    | Shock sensor threshold                 | YES |
|              | R                | Read Shock sensor threshold from panel | NO  |
|              | W                | Write Shock sensor threshold to panel  | NO  |
| Syncbox      | Tube Ready Time  | Not used                               | NO  |
|              | R                | Not used                               | NO  |
|              | W                | Not used                               | NO  |
| Button       | L                | Read Log from panel                    | NO  |


#### 4.13.3.3 Network Settings

Client Mode:



AP Mode:

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 107/137      |




| Item        |             | Description  | Modify |
|-------------|-------------|--|--------|
| Mode        | Client      | Set panel in client mode                                     | NO     |
|             | AP          | Set panel in AP mode   | NO     |
| Client Mode | Add         | Add available wireless AP account                            | NO     |
|             | Del         | Delete Exist wireless AP account                             | NO     |
|             | Up          | Wireless AP account move up                                  | NO     |
|             | Down        | Wireless AP account move down                                | NO     |
|             | Default     | Set AP account as default connection                         | NO     |
|             |             |  |        |
| AP Mode     | SSID        | Wireless AP SSID when panel in AP mode                       | YES    |
|             | Key         | Wireless AP Key when panel in AP mode                        | YES    |
|             | Security    | Wireless AP Security method when panel in AP mode            | YES    |
|             | Frequency   | Wireless AP frequency(2.4GHz and 5GHz) when panel in AP mode | YES    |
|             | Country     | Wireless AP Country when panel in AP mode                    | YES    |
|             | Band        | Wireless AP Band(HT20 and HT40) when panel in AP mode        | YES    |
|             | Channel     | Wireless AP Channel when panel in AP mode                    | YES    |
|             | DHCP Server | DHCP function when panel in AP mode                          | YES    |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 108/137      |

|        |                      |  |     |
|--------|----------------------|--|-----|
|        | Scan Local Wlan      | Scan local wifi signal when panel in AP mode | YES |
| FPD    | IP                   | Network IP address of panel                  | YES |
|        | Port                 | Network Port of panel                        | NO  |
|        | MAC                  | Network MAC address of panel                 | YES |
|        | Image Channel        | Network protocol of panel                    | YES |
| Button | Read Wifi Status     | Read wireless module status from panel       | NO  |
|        | Read Config from FPD | Read wireless module setting from panel      | NO  |
|        | Write Config to FPD  | Write wireless module setting to panel       | NO  |
|        | Scan from FPD        | Scan Wifi signal by panel                    | NO  |

#### 4.13.3.4 HW Correction




| Item | Description   | Modify |
|------|---|--------|
| HW   | Offset Mode<br>Hardware offset mode of panel<br>NO: no hardware offset mode<br>Pre: hardware pre-offset mode<br>Post: hardware post-offset mode | YES    |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 109/137      |

|  |               |  |    |
|--|---------------|--|----|
|  | Gain Map      | Gain calibration mode of panel<br>Do in SDK: software gain calibration<br>Do in HW: hardware gain calibration    | NO |
|  | Defect Map    | Defect correction mode of panel<br>Do in SDK: software defect correction<br>Do in HW: hardware defect correction | NO |
|  | Lag Map       | Lag correction of panel  | NO |
|  | Download      | Download correction and calibration template to panel  | NO |
|  | Read out      | Upload correction and calibration template from panel  | NO |
|  | Update local  | Replace local correction and calibration template with template uploaded currently                               | NO |
|  | Update status | Get correction and calibration configuration from panel  | NO |

#### 4.13.3.5 Firmware Update



Config

**Firmware Update**

Status:

Firmware

Update Type:

Image File:

Progress:

Recover ROM

Write ROM

Read ROM

Write RAM

Read RAM

Reset FPD

Exit

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 110/137      |

|          | Item        | Description                  | Modify |
|----------|-------------|------------------------------|--------|
| Firmware | Update Type | Not used                     | YES    |
|          | Image File  | Local address of update file | NO     |
|          | Progress    | Progress bar of updating     | NO     |
|          | Select      | Select update file           | NO     |
|          | Update      | Start update                 | NO     |
|          | Reload      | Reload firmware of panel     | NO     |

## 4.13.4 Correction and Calibration

### 4.13.4.1 Gain Calibration

**Correction Wizard**

Gain | Defect

PREP   Acquire   PREP Acquire

Expected Average

Expected Gray Value

5500

=>

New Point

Selected Gain 0

Current Average 0

Center Average 0

=>

Selected Points

New Gain Point

Select

Create Gain Correct File


OK

|                  | Item                | Description                         | Modify |
|------------------|---------------------|-------------------------------------|--------|
| Expected Average | Expected Gray Value | Not used                            | YES    |
|                  | New Gain Point      | Start generating gain template      | NO     |
| New Point        | Selected Gain       | Number of selected gain point       | NO     |
|                  | Current Average     | Average value of the complete image | NO     |

|  |  |         |            |
|--|--|---------|------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02 |
|  |  | Version | A2         |
|  | Mars1717V User' s Manual               | Date    | 2016.02.23 |
|  |  | Page    | 111/137    |

|        |                          |  |    |
|--------|--------------------------|--|----|
| Button | Center Average           | Average value of the image in center ROI<br>( 100X100) | NO |
|        | Select                   | Select and save current gain point                     | NO |
|        | Prep                     | Send “Clear” command to panel                          | NO |
|        | Acquire                  | Send “Acquire” command to panel                        | NO |
|        | Prep Acquire             | Send “Clear Acquire” command to panel                  | NO |
|        | Create Gain Correct File | Start generating gain Calibration template             | NO |
|        |                          |  |    |
|        |                          |  |    |

#### 4.13.4.2 Defect Correction




| Item              |                | Description                      | Modify |
|-------------------|----------------|----------------------------------|--------|
| Create Light File | Current        | Current sequence of defect point | NO     |
|                   | Start creating | Start defect point acquisition   | NO     |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 112/137      |

|        |                   |                                       |    |
|--------|-------------------|---------------------------------------|----|
|        | Auto completion   | Automatically acquire defect point    | NO |
| Button | Refresh HW Offset | Updating pre-offset template in panel | NO |
|        | Template          |                                       |    |
|        | Prep              | Send “Clear” command to panel         | NO |
|        | Acquire           | Send “Acquire” command to panel       | NO |
|        | Prep Acquire      | Send “Clear Acquire” command to panel | NO |
|        | Create Offset     | Generate pre-offset template          | NO |
|        | Create Defect     | Generate Defect correction template   | NO |

#### 4.13.5 Image Check and upload




|                  | Item        | Description                                   | Modify |
|------------------|-------------|---|--------|
| Remote File List | Index       | Sequence number of image                      | NO     |
|                  | File name   | Name of image in panel                        | NO     |
|                  | Create Time | Acquisition time of image in panel            | NO     |
|                  | Attr        | Image label<br>0x01—Do Offset<br>0x02—Do Gain | NO     |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 113/137      |

|        |                 |   |    |
|--------|-----------------|---|----|
|        |                 | 0x04—Do Most gain<br>0x08—Do Defect<br>0x10—Post Offset Raw Image |    |
|        | Delay Time      | Time between clear and image acquisition                          | NO |
| Button | Get File List   | Get file list from panel  | NO |
|        | Start Upload    | Start uploading image selected                                    | NO |
|        | Stop Upload     | Stop uploading process  | NO |
|        | Open in List    | Open image selected   | NO |
|        | Open Local File | Open local images in Workstation                                  | NO |

#### 4.13.6 Infrared Registration



| Item | Description |  | Modify |
|------|-------------|--|--------|
| /    | Work Mode   | Work mode of infrared registration tools<br>Write Only: infrared registration tools is allowed to write to panel | YES    |



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 114/137      |

|                           |           |   |     |
|---------------------------|-----------|---|-----|
|                           |           | <p>Read Only: infrared registration tools is allowed to read from panel</p> <p>Read &amp; Write: infrared registration tools is allowed to read from panel and write to panel</p> <p>Read &amp; confirm by User: infrared registration tools is allowed to read from panel and write to panel only when confirmed by user</p> |     |
| AP Mode Configuration     | AP Mode   | Set panel in AP mode or Client mode   | YES |
|                           | SSID      | Wireless AP SSID when panel in AP mode  | YES |
|                           | Key       | Wireless AP Key when panel in AP mode   | YES |
|                           | Security  | Wireless AP Security way when panel in AP mode  | YES |
|                           | Frequency | Wireless AP Frequency(2.4GHz and 5GHz) when panel in AP mode  | YES |
|                           | Country   | Wireless AP Country Code when panel in AP mode  | YES |
|                           | Band      | Wireless AP Band(HT20 and HT40) when panel in AP mode   | YES |
|                           | Channel   | Wireless AP Channel when panel in AP mode   | YES |
| Client Mode Configuration | SSID      | Wireless SSID when panel in Client mode   | YES |
|                           | Key       | Wireless Key when panel in Client mode  | YES |
| Button                    | Apply     | Save wireless parameter in infrared registration tools  | NO  |
|                           | Start     | Start write wireless parameter in panel   | NO  |
|                           | Exit      | Exit infrared registration tools  | NO  |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 115/137      |

## 5 Regulatory Information

### 5.1 Medical equipment safety standards

#### ◆ Medical equipment classification

|   |  |
|---|--|
| Type of protection against electrical shock   | Class I Equipment, using medical approved adaptor supply<br>Internally powered Equipment, using battery power supply   |
| Degree of protection against electrical shock | Without Applied Parts  |
| Degree of protection against ingress of water | IPX0 (Mars1717V)<br>IPX0 (Charger-KV)  |
| Mode of operation                             | Continuous operation   |
| Flammable anesthetics                         | Not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide<br>Not suitable for use in the oxygen rich environment |

Note: The product safety standards that apply to Mars1717V which includes the main units: detector and charger-kv .

#### ◆ References harmonized standards under Directive 93/42/EEC

|   |  |
|---|--|
| MDD (93/42/EEC)                             | Medical Device Directive   |
| EN ISO 13485:2012/EN ISO 13485:2012/AC:2012 | Medical devices --- Quality management systems --- Requirements for regulatory purposes                              |
| EN ISO14971: 2012                           | Medical device – Application of risk management to medical devices   |
| EN 60601-1:2014                             | Medical electrical equipment -- Part 1: General requirements for basic safety and essential performance              |
| EN 60601-1-2:2007                           | Medical electrical equipment – Part 1-2: Collateral standard: Electromagnetic compatibility – Requirements and tests |
| EN 60601-2-54:2009                          | Medical electrical equipment -- Part 2-54: Particular requirements for the basic safety and essential performance of |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 116/137      |

|                       |  |
|-----------------------|--|
|                       | X-ray equipment for radiography and radioscopy   |
| IEC 62133:2012        | Secondary cells and batteries containing alkaline or other non-acid electrolytes – Safety requirements for portable sealed secondary cells, and for batteries made from them, for use in portable applications |
| EN 62220-1:2004       | Medical electrical equipment - Characteristics of digital X-ray imaging devices - Part 1: Determination of the detective quantum efficiency  |
| EN 62304:2006/AC:2008 | Medical device software - Software life-cycle processes  |
| EN 62366:2008         | Medical devices - Application of usability engineering to medical devices  |

## 5.2 Guidance and manufacture's declaration for EMC

### ◆ Electromagnetic emissions

Mars1717V is intended for use in the electromagnetic environment specified below. The customer or the user of Mars1717V should assure that it is used in such an environment.

| Emission Test   | Compliance | Electromagnetic Environment - Guidance   |
|---|------------|--|
| RF emissions CISPR 11                                       | GROUP1     | Mars1717V uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.<br><br>Mars1717V is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| RF emissions CISPR 11                                       | Class B    |  |
| Harmonic emissions<br>IEC 61000-3-2                         | Class B    |  |
| Voltage fluctuations/ flicker<br>emissions<br>IEC 61000-3-3 | Complies   |  |

### ◆ Electromagnetic immunity

|  |  |                |              |
|--|--|----------------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | <b>Doc NO.</b> | 037-201-02   |
|  |  | <b>Version</b> | A2           |
|  | Mars1717V User' s Manual               | <b>Date</b>    | 2016. 02. 23 |
|  |  | <b>Page</b>    | 117/137      |

Mars1717V is intended for use in the electromagnetic environment specified below. The customer or the user of Mars1717V should assure that it is used in such an environment.


| Immunity Test  | IEC 60601 Test Level  | Compliance Level  | Electromagnetic Environment - Guidance  |
|--|---|---|---|
| Electrostatic discharge (ESD)<br>IEC 61000-4-2   | ±6 kV contact<br>±8 kV air  | ±6 kV contact<br>±8kV air   | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.   |
| Electrical fast transient/ burst<br>IEC 61000-4-4  | ±2 kV for power supply lines<br>±1 kV for input/output lines  | ±2 kV for power supply lines<br>±1 kV for input/output lines  | Mains power quality should be that of a typical commercial or hospital environment.   |
| Surge<br>IEC 61000-4-5   | ±1 kV differential mode<br>±2 kV common mode  | ±1 kV differential mode<br>±2 kV common mode  | Mains power quality should be that of a typical commercial or hospital environment.   |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC 61000-4-11 | <5% UT (>95% dip in UT) for 0.5 cycle.<br>40% UT (60% dip in UT) for 5 cycle.<br>70% UT (30% dip in UT) for 25 cycle.<br><5% UT (>95% dip in UT) for 5 sec. | <5% UT (>95% dip in UT) for 0.5 cycle.<br>40% UT (60% dip in UT) for 5 cycle.<br>70% UT (30% dip in UT) for 25 cycle.<br><5% UT (>95% dip in UT) for 5 sec. | Mains power quality should be that of a typical commercial or hospital environment. If the user of that requires continued operation during power supply interruptions, it is recommended that be powered from an uninterruptible power supply. |
| Power frequency (50/60Hz) magnetic field<br>IEC 61000-4-8  | 3 A/m   | 3 A/m   | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.   |

UT is the a.c. mains voltage prior to application of the test level.

#### ◆ Guidance and manufacturer's declaration---electromagnetic immunity

Mars1717V is intended for use in the electromagnetic environment specified below. The customer or the user of Mars1717V should assure that it is used in such an environment.

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 118/137      |

| Immunity Test                 | IEC 60601 Test Level      | Compliance Level | Electromagnetic Environment - Guidance  |
|-------------------------------|---------------------------|------------------|---|
| Conducted RF<br>IEC 61000-4-6 | 3 Vrms<br>150kHz to 80MHz | 3 Vrms           | Portable and mobile RF communications equipment, AC-DC adapter or electromagnet should be used not closer to any part of the Model Mars1717V, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.<br><br>Recommended separation distance<br><br>$d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 2,3 \sqrt{P} \text{ 800 MHz to 2,5 GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).<br><br>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.<br><br>Interference may occur in the vicinity of equipment marked with the following symbol:  |
| Radiated RF<br>IEC 61000-4-3  | 3 V/m<br>80MHz to 2.5GHz  | 3 V/m            |   |

NOTE: UT is the a.c. mains voltage prior to application of the test level.

NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Mars1717V is used exceeds the applicable RF compliance level above, Mars1717V should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating Mars1717V.

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 119/137      |

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

- ◆ Recommended separation distances between portable or mobile RF communications equipment and Mars1717V

Mars1717V is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of Mars1717V can help prevent electromagnetic interference by maintaining a minimum distance between portable or mobile RF communications equipment (transmitters) and Mars1717V as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter /W | Separation distance according to frequency of transmitter /m |                                      |                                       |
|--|--|--------------------------------------|---------------------------------------|
|  | 150kHz~80 MHz<br>$d = 1.2 \sqrt{P}$                          | 80 MHz~800 MHz<br>$d = 1.2 \sqrt{P}$ | 800 MHz ~2.5GHz<br>$d = 2.3 \sqrt{P}$ |
| 0.01   | 0.12   | 0.12                                 | 0.23                                  |
| 0.1  | 0.38   | 0.38                                 | 0.73                                  |
| 1  | 1.2  | 1.2                                  | 2.3                                   |
| 10   | 3.8  | 3.8                                  | 7.3                                   |
| 100  | 12   | 12                                   | 23                                    |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- ◆ Cables information below is provided for EMC reference.

| Cable | Recommended cable length | Shielded or Unshielded | Number | Cable classification |
|-------|--------------------------|------------------------|--------|----------------------|
|-------|--------------------------|------------------------|--------|----------------------|

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 120/137      |

| Cable             | Recommended cable length | Shielded or Unshielded | Number | Cable classification |
|-------------------|--------------------------|------------------------|--------|----------------------|
| Input Power Cable | 3m                       | Unshielded             | 1 Set  | AC Power             |
| Ethernet Cable    | 3.5m                     | Shielded               | 1 pcs  | Signal               |
| DC Power Cable    | 3.5m                     | Unshielded             | 1 pcs  | Signal               |
| LAN Cable         | 3m                       | Shielded               | 1 pcs  | Signal               |

◆ Important information regarding Electro Magnetic Compatibility (EMC)

Mars1717V require special precautions regarding EMC and needs to be installed only by iRay or authorized personnel and put into service according to EMC information provided in the user manual. Mars1717V in use may be susceptible to electromagnetic interference from portable and mobile RF communications such as mobile (cellular) telephones. Electromagnetic interference may result in incorrect operation of the system and create a potentially unsafe situation.

Mars1717V conforms to this EN60601-1-2:2007 standard for both immunity and emissions.

Nevertheless, special precautions need to be observed:

The use of accessories, transmitters and cables other than those specified by this User Manual, with the exception of accessories and cables sold by iRay of Mars1717V as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of Mars1717V.

Mars1717V should not be used adjacent to or stacked with other equipment.

### 5.3 Radio Frequency Compliance Information

| Country | Item  |
|---------|---|
| U.S.A   | FCC Part 15.107 Subpart (b) / 15.109(g) Subpart B<br>FCC Part 15 Subpart E 15.407<br>FCC Part 15 Subpart C 15.247 |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 121/137      |

|                |  |
|----------------|--|
| European Union | ETSI EN 301 489-1 V1.8.1 (EMC)<br>ETSI EN 301 489-17 V2.1.1 (EMC)<br>EN 300 328 V1.7.1; EN 301 893 V1.6.1 (RF)<br>EN 62311:2008 (RF Exposure)<br>ETSI EN 300 328 V1.7.1; EN 301 893, V1.5.1 (Radio Spectrum) |
|----------------|--|

### 5.3.1 FCC Compliance

- This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.
- Operation is subject to the following tow conditions.

This device may not cause harmful interference.

This device must accept any interference received, including interference that may cause undesired operation.

- This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measure.

Reorient or relocate the receiving antenna.

Increase the separation between the equipment and receiver.

Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

Consult the distributor or an experienced radio/TV technician for help.

## 5.4 Battery Safety Standards

| Standards | Description  |
|-----------|--|
| UL1642    | Component Recognition on the Secondary Li-ion cell |




|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 122/137      |

|                    |   |
|--------------------|---|
| UL 2054:2004 R9.11 | Household and commercial Batteries  |
| IEC 62133:2012     | Secondary cells and batteries containing alkaline or other non-acid electrolytes  |
| UN38.3             | United Nations Recommendations on the Transport of dangerous goods Manual of tests and Criteria<br>ST/SG/AC.10/11/Rev.5/Amend.1&Amend.2 |

## 5.5 Product Label

### 5.5.1 Detector

Mars1717V-PSI detector



Mars1717V-VSI detector

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 123/137      |

## 无线数字平板探测器

产品型号: Mars1717V-VSI  
接入电源: 适配器供电输入 24V  $\equiv$  0.8A  
电池供电输入 10.8V  $\equiv$  1.78A  
附属设备 Mars1717V-VSI GB9706.14-1997

上海奕瑞光电科技有限公司  
中国上海浦东张江东区医疗器械园  
瑞庆路590号9幢2层202室

医疗器械注册证编号: 沪械注准201X231XXXX  
其它内容详见说明书



Model: Mars1717V-VSI  
Power: Adapter Port Input 24V  $\equiv$  0.8A  
Battery Port Input 10.8V  $\equiv$  1.78A  
iRay Technology ( Shanghai ) Ltd.  
RM202, Building 7, No. 590, Ruiqing RD.,  
Zhangjiang East, Pudong, Shanghai, China  
www.iraychina.com

EC REP Renault - Petersen Limited  
COUCHING HOUSE, COUCHINE STREET,  
WATLINGTON, OXFORDSHIRE, OX49 5PX, UK

FCC ID: 2ACHK-02110113



CE REP

FC

SN



\*\*\*\*\*


CE 0197




|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 124/137      |

## 5.5.2 Battery

### 电池

产品型号: BATTERY-KV  
 输出: 10.8V === 4180mAh/45.14Wh  
 充电电压: 12. 6V  
 执行标准: GB/T18287-2013

 上海奕瑞光电子科技有限公司  
 中国上海浦东张江东区医疗器械园  
 201X-XX 瑞庆路590号9幢2层202室


**警告!**


- 注意火灾、爆炸危险，禁止灼烧或浸入水中；
- 不可短接、不可存放于高于45℃的环境或者擅自拆解；
- 请使用CHARGER-KV为电池充电，请勿使用其他型号充电器；
- 请在初次使用时请为电池完全充电；
- 当电池电量小于5%时，请及时给电池充电；
- 请用同型号电池替代，使用不同型号电池有燃烧或者爆炸危险；
- 请将电池放置在儿童无法触及的地方；
- 请根据当地法律处理废旧电池。

---

### Battery





Model: BATTERY-KV  
 Output: 10.8V === 4180mAh/45.14Wh  
 Charge Voltage: 12. 6V  
 Standard: GB/T18287-2013

 iRay Technology (Shanghai) Ltd.,  
 RM202, Building 7, No.590, Ruiqing RD.,  
 Zhangjiang East, Pudong, Shanghai, China  
 201X-XX


 Renault - Petersen Limited  
 COUCHING HOUSE, COUCHINE STREET,  
 WATLINGTON, OXFORDSHIRE, OX49 5PX, UK

**CAUTION:**

- Risk of fire, explosion, burning or put into water.
- Do not short circuit, crush or expose battery in the environment above 45℃, incinerate or disassemble the battery.
- Charge only with CHARGER-KV.
- Please charge the battery before first use.
- Please charge the battery immediately when the capacity is less than 5%.
- Replace battery with same part number only. Use another battery may present a risk of fire or explosion.
- Keep away from children.
- Dispose of all used batteries according to local law.

SN



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 125/137      |

### 5.5.3 Battery Charger



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 126/137      |

## 6 Trouble Shooting

When user encounters problems or error messages, refer to this chapter. If the problem persists, turn off the panel and contact iRay service department ([service@iraychina.com](mailto:service@iraychina.com)). We would provide the best service.

### 6.1 Cable Inspection

#### 6.1.1 Ethernet cable


##### 6.1.1.1 Inspection method

Test the Ethernet cable detector by cable tester, and confirm whether all the cores of the cable are conductive. If the cable tester is not available, please check the definition of Ethernet cable and test the conductive between each pin of connector.




##### 6.1.1.2 Definition of Ethernet cable

| Colour         | J1    | J2    |
|----------------|-------|-------|
| Orange\White   | 1     | 3     |
| Orange         | 2     | 6     |
| Green\White    | 3     | 1     |
| Blue           | 4     | 7     |
| Blue\White     | 5     | 8     |
| Green          | 6     | 2     |
| Brown\White    | 7     | 4     |
| Brown          | 8     | 5     |
| Grounding wire | Shell | Shell |



##### 6.1.1.3 Charge cable inspection


| J1 | J2 |
|----|----|
| 3  | 1  |
|    | 3  |
| 2  | 2  |
|    | 4  |



Measure the conductive between each pin of J1 and J2 by multimeter.

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 127/137      |


### 6.1.2 Date cable inspection



According to the definition of data cable to verify the conductive between each pin of J1 and J2 by multimeter if there is only one data cable in the field.


## 6.2 Adapter Inspection

### 6.2.1 Adapter connector definition



| Pin | Description  |
|-----|--------------|
| 2   | DC +24V/5.0A |
| 3   | 0V           |


## 6.2.2 Inspection method



Connect the adapter to the detector or charging dock to confirm whether the battery is in working order if verifying the output current is not practicable.

## 6.3 Battery Inspection


### 6.3.1 Battery pin definition



| Pin | Symbol | Description                         |
|-----|--------|-------------------------------------|
| 1,2 | P+     | Battery Discharge Positive Terminal |
| 7,8 | P-     | Battery Discharge Negative Terminal |

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 129/137      |


## 6.3.2 Inspect Method



The battery is expendable and its life is shortened through use. Recycle the battery if the fully charged battery only last not more than 1 hour. Wasted batteries suggested return to manufacturer or put at appointed public battery reclaim area, do not mix battery with other waste or dispose of battery ad libitum.

## 6.4 Dock Inspection


### 6.4.1 Dock pin definition





| Pin | Voltage | Description                      |
|-----|---------|----------------------------------|
| 1,2 | +       | Battery Charge Positive Terminal |
| 7,8 | -       | Battery Charge Negative Terminal |

## 6.4.2 Inspection method




## 6.5 Detector Main Unit Inspection

Know how to obtain the log of SDK and detector is necessary before diagnose any problem of detector main unit, and the method is described in chapter 6.1.

### 6.5.1 Get SDK and detector log

#### 6.5.1.1 Fetch SDK log

Find the location of iDemo.exe, and there is a folder called “Correct” in the same directory. Find the folder named by the serial number of detector. The folder called “Log” is the storage path of SDK fie, and the date is contained in name of each log file. Please compress the log file before send it.





The log function can be switched on in the config.ini. in the folder named “Res”, please check the configuration before start iDemo.


SET\_LOG\_INFO=1//Open log function'  
SET\_LOG\_INFO=0//Close log function

```
// *** log *** //
SET LOG INFO=1
```


|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 131/137      |

### 6.5.1.2 Detector log


Find FTP Server.exe  and open the FTP Server, attention that the FTP port number should be set as same as the config.ini of iDemo. Click  button to start FTP server after finish setting the FTP port number.



Find the “L” button in the following diagram in the configuration GUI of iDemo, click “L” button and the FTP will upload the “Log” in the storage path “upload\[FPD 13 number S/N]\logs” which is at the same location of iDemo.exe.



## 6.6 Power up failure inspection




is the label of power indicator. Skip the step “Press reset button” if no reset button, because the previous version doesn’t content press button, and it’s no effect in normal use.

## 6.7 Connection inspection




is label of link indicator. If there is connection error from iRayDR, please view the FAQ of iRay before checking the connection issue through iDemo

### 6.7.1 Connection failed in wired mode







## 6.7.2 Connection Failed in Client Mode



## 6.7.3 Disconnect in AP Mode




## 6.8 Dark Image Acquisition After Exposure in Inner Mode



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 134/137      |

## 6.9 No Image Acquire after Exposure



## 6.10 Image Inspection


All the correction files are in the storage path “res\Correct\ [FPD series number] \0”

| Calibration type | File format | Remark   |
|------------------|-------------|--|
| Offset           | *.off       | iRay mostly use post offset, the offset is done during the image acquisition process. So the offset file is not effective in the correction file storage path. |
| Gain             | *.gn        | The gain file can be select or deselect through iDemo  |
| Defect           | *.dft       | The gain file can be select or deselect through iDemo  |
| Lag              | *.lag       | The Most gain is effective while it's in the correction file storage path. Rename it or move it to other directory can make it invalid.                        |

The following is the image

| Image type               | Description  |
|--------------------------|--|
| Original dark image      | Acquire by click “Prep Acquire” button without exposure and offset correction                                |
| Dark image               | Acquire by click “Prep Acquire” button without exposure and the offset calibration should be added at least. |
| Original bright image    | Acquire the image under exposure condition and without any correction.                                       |
| Bright image with offset | Acquire the image under exposure condition and add offset correction only.                                   |


### 6.10.1 Artifact on bright image




### 6.10.2 Horizontal Line on bright image

Adjust the WW to 30 and WL to 100 through iDemo for checking the horizontal line of dark image. Normally the dark image with horizontal line is shown as the figure 1.

|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 135/137      |




### Inspection method




|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A2           |
|  | Mars1717V User' s Manual               | Date    | 2016. 02. 23 |
|  |  | Page    | 136/137      |

### 6.10.3 Defective line




### 6.10.4 Defective point



### 6.10.5 TFT broken

Acquire a dark image, if one part of the image is much smooth than the normal dark image. The TFT must be broken by intense impact.





|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A0           |
|  | Mars1717V User's Manual                | Date    | 2015. 08. 26 |
|  |  | Page    | 137/137      |

## 7 Service Information

### 7.1 Product lifetime

The estimated product lifetime is up to 10 years under appropriate regular inspection and maintenance.

### 7.2 Regular inspection and Maintenance

In order to ensure the safety of patients, operating person and third parties, and to maintain the performance and reliability of the equipment, be sure to perform regular inspection at least once a year. If necessary, clean up the equipment, make adjustments, or replace consumables such as fuses, detector cable, etc. There may be cases where overhaul is recommended depending on conditions. Contact iRay service office or local iRay dealer for regular inspection or maintenance.

There is a Ni-MH battery in the FPD, its lifetime is 5 years, when arrived in the lifetime of the battery is need to be placed. And the placement need contact Shanghai Iray after-sales service departments or authorized product distributors.



|  |  |         |              |
|--|--|---------|--------------|
|  | <b>iRay Technology (Shanghai) Ltd.</b> | Doc NO. | 037-201-02   |
|  |  | Version | A0           |
|  | Mars1717V User's Manual                | Date    | 2015. 08. 26 |
|  |  | Page    | 138/137      |

## 7.3 Repair

If a problem cannot be solved even taking the measures indicated in troubleshooting, contact your sales representative or local iRay dealer for repairs. Please refer to the name label and provide the following information:

Product Name:

Series Number:

Description of Problem: as clearly as possible.

## 7.4 Replacement parts support

Performance parts (parts required to maintain the function of the product) of this product will be stocked for 5 years after discontinuance of production, to allow for repair.